# **LifePort<sup>®</sup>** **INFUSION SET**

INSTRUCTIONS FOR USE

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MODE D'EMPLOI

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**ANGIODYNAMICS<sup>®</sup>**

# ANGIODYNAMICS®

## TABLE OF CONTENTS


<b>en</b>	<b>INSTRUCTIONS FOR USE</b>	<b>1</b>
<b>du</b>	<b>GEBRUIKSAANWIJZING</b>	<b>8</b>
<b>es</b>	<b>INSTRUCCIONES DE USO</b>	<b>15</b>
<b>it</b>	<b>ISTRUZIONI PER L'USO</b>	<b>23</b>
<b>fr</b>	<b>MODE D'EMPLOI</b>	<b>30</b>
<b>gm</b>	<b>GEBRAUCHSANWEISUNG</b>	<b>38</b>

# **LifePort**<sup>®</sup>

## **Infusion Set**

**CAUTION:** Federal (U.S.A.) law restricts this device to sale by, **RX only** or on the order of, a physician.

**Figure 1**

**ATTENTION:** See Instructions for Use

STERILE EO

Inner package contents sterile, if not opened or damaged. Sterilized using Ethylene Oxide

Non-Pyrogenic. Do Not Reuse. Single Use Only

## INDICATIONS FOR USE:

The LifePort® Infusion Set is used to access implanted vascular ports to administer fluids and/or to withdraw blood.

## MRI INFORMATION:

The LifePort® Infusion Set was determined to be MR-conditional according to the terminology specified in the American Society for Testing and Materials (ASTM) International, Designation: F2503-05. Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment. ASTM International, 100 Barr Harbor Drive, PO Box C700, West Conshohocken, Pennsylvania, 2005.

Non-clinical testing demonstrated that the LifePort® Infusion Set is MR Conditional. A patient with this device can be scanned safely immediately after placement under the following conditions:

- Static magnetic field of 3-Tesla or less
- Spatial gradient magnetic field of 720-Gauss/cm or less

## MRI-Related Heating

In non-clinical testing, the LifePort® Infusion Set produced the following temperature rise during MRI performed for 15-min in the 3-Tesla (3-Tesla/128-MHz, Excite, Software G3.0-052B, General Electric Healthcare, Milwaukee, WI) MR system:

*Highest temperature change +1.7° C*

Therefore, the MRI-related heating experiments for the LifePort® Infusion Set at 3-Tesla using a transmit/receive RF

body coil at an MR system reported whole body averaged SAR of 3.0-W/kg (i.e., associated with a calorimetry measured value of 2.8-W/kg) indicated that the greatest amount of heating that occurred in association with these specific conditions was equal to or less than +1.7° C.

### Artifact Information

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the LifePort® Infusion Set. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary.

Pulse Sequence	T1-SE	T1-SE	GRE	GRE
Signal Void Size	1,299-mm <sup>2</sup>	1,130-mm <sup>2</sup>	4,987-mm <sup>2</sup>	4,081-mm <sup>2</sup>
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular

### CONTRAINDICATIONS:

Do not use with patients showing signs or symptoms of infection at the insertion site.

Do not use with patients who have exhibited prior intolerance to the materials of construction.

### WARNINGS:

Maximum pressure recommended is 40 psi. A 10 ml syringe or larger is recommended.

Do not use the device if the package has been opened or damaged.

Clamp infusion set when not in use.

Store in a cool dry place.

**CAUTIONS:**

**RX only** Federal (USA) law restricts this device to sale by or on the order of a physician.

**SINGLE USE ONLY:** Non-Pyrogenic, Non-toxic.  
Contents sterile unless package is opened or damaged.  
Do not resterilize.

See instructions for use.

Carefully read and follow all instructions prior to use of the device.

Follow CDC guidelines, institutional guidelines or INS and ONS guidelines for infusion set use.

Always maintain universal precautions and utilize aseptic technique throughout insertion care and maintenance procedures.

Application of excessive force to needle while accessing port may result in damage to needle and port.

**KEEP HAND / FINGERS PROTECTED FROM NEEDLE TIP AT ALL TIMES BY FOLLOWING THE INSTRUCTIONS FOR USE.**

**INSTRUCTIONS:****PROCEDURE FOR PORT ACCESS:**

**NOTE:** Aseptic technique should be maintained during port access.

**NOTE:** Follow Institutional Universal Precautions.

1. Gather supplies necessary for port access.
2. Place mask on face.
3. Open and drop all supplies onto sterile field.
4. Don sterile gloves.
5. Cleanse area around port, following your institutional protocol.
6. Remove luer cap and attach injection cap used by your facility.

**NOTE:** Prior to infusing, ensure that injection cap is securely connected.

7. Prime infusion set following your institutional protocol.

## **INSERTION:**

1. Position index finger at the bend of the needle.
2. Remove protective sleeve on needle.
3. Insert needle at a 90° angle to the port.

**NOTE:** For optimum needle stability during insertion, hold wings and infusion set body as illustrated in Figure 1.

4. Advance needle through the skin and the port septum until it contacts the bottom of the port reservoir.  
Excessive force may lead to bending of the bevel and cause pain and difficulty with removal of the needle from the septum.

**NOTE:** Needle placement should be confirmed by aspiration of 3-5 ml of blood. If blood return is not present, the implanted port should be evaluated by the physician.

**NOTE:** A 10 ml syringe or larger is recommended.

5. Secure the device per institutional protocol.
6. Using turbulent “push-pause” technique, flush port per institutional protocol maintaining positive pressure.

### **REMOVAL:**

1. Using turbulent “push-pause” technique, flush the non-coring needle per institutional protocol, maintaining positive pressure.
2. Remove dressing.
3. Remove needle. Use caution to avoid rebound injury.
4. Properly dispose in Sharps Container.

### **INJECTION CAP CHANGE RECOMMENDATION:**

Change injection caps per institutional guidelines or when cap has been removed for any reason, appears damaged, is leaking, or if there is blood residue on the cap.

### **INJECTION CAP CHANGE PROCEDURE:**

1. Wash hands thoroughly with soap and warm water.
2. Don clean gloves.
3. Ensure tubing clamp is closed prior to disconnecting the injection cap.
4. Open the sterile injection cap package.
5. Purge air from device by flushing with appropriate solution.
6. Clean the tubing’s luer lock hub on the outside with povidone-iodine or swab with alcohol. Allow to dry.
7. Remove the protective piece from the new injection cap, being careful not to touch the inner portion of the new injection cap or infusion set luer.
8. Attach the new injection cap to the luer in a clockwise manner.  
**DO NOT OVER TIGHTEN.**

- \* Approximate priming volumes:  
Infusion set with y-site: 0.6 ml  
Infusion set without y-site: 0.4 ml

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LifePort is a registered trademark of  
AngioDynamics, Inc.  
under U.S. Pat# 4769010.

SEE INSTRUCTIONS FOR USE

SINGLE USE ONLY, DO NOT REUSE. DO  
NOT RESTERILIZE

STERILE EO

STERILIZED WITH ETHYLENE OXIDE.  
NON-PYROGENIC. DO NOT RESTERILIZE.  
STERILE IF PACKAGE IS NOT OPENED OR DAMAGED.

CAUTION: Federal (U.S.A.) Law restricts this  
device to sale by or on the order of a physician.

ANGIODYNAMICS®

**Manufactured by:**  
**AngioDynamics, Inc.**  
One Horizon Way  
Manchester, GA 31816  
P: 518.798.1215 or  
800.772.6446

**Authorized Representative:**  
**AngioDynamics UK Ltd.**  
Building 2000, Beach Drive  
IQ Cambridge, Waterbeach,  
Cambridge, CB25 9TE  
United Kingdom  
Telephone: +44 (0) 122-372-9361  
Fax: +44 (0) 122-372-9329

[www.angiodynamics.com](http://www.angiodynamics.com)