**EN INDICATIONS**  
**ANGIODYNAMICS** Angiographic Catheters are for use where angiographic diagnosis is indicated.

#### POTENTIAL ADVERSE EFFECTS

The following adverse reactions have been reported and are associated with the use of angiographic catheters:

- Thrombus formation
- Emboli
- Arterial wall damage
- Plaque dislodgment
- Hematoma
- Cardiac arrhythmias
- Myocardial infarction
- Stroke
- Death

#### WARNINGS

- ANGIODYNAMICS Angiographic Catheters should be used only by physicians with a thorough understanding of angiography and percutaneous interventional procedures and are not for reuse or sterilization.
- ANGIODYNAMICS 3 French Angiographic Catheters are not intended for intracranial use.
- Do not insert catheters directly through synthetic vascular grafts. Insert through a sheath introducer.
- Never advance or retract an Angiographic Catheter or guidewire against resistance. This may result in damage to the vessel, the product, or both.
- Do not attempt to hand straighten the tip of any curved tipped catheters which are furnished with a tip straightener. This may result in damage to the product. Tip straighteners are furnished on all catheters intended to be straightened with the aid of a tip straightener.
- Reshaping of the catheter tip is not recommended. Physical damage to the catheter material can result when exposed to heat.
- Do not leave curved tipped catheters straightened over guidewires for extended periods of time. This will result in failure of the catheter tip to re-form to its intended shape.
- When using ANGIODYNAMICS HALO Catheters, form or straighten the tip in the distal thoracic aorta using a guidewire when advancing or withdrawing the catheter.
- Avoid injection of contrast medium into the intercostal and lumbar arteries.
- Always use a guidewire to remove the catheter from the vasculature. Failure to do so may result in damage to the vessel, puncture site, product, or all three.
- When retracting catheters, great care must be taken to avoid exerting excessive pressure at the groin entry site. Excessive pressure may result in damage to the vessel, the catheter, or both.
- The maximum pressure limit of catheters intended for flush angiography is stated on the catheter package. When using a pressure injector, do not exceed the stated maximum pressure.
- Catheters intended for selective angiography do not have a maximum pressure limit stated on the catheter package. Typical flow rates of up to 10cc per second are stated with pressure generated to achieve these typical flow rates.
- Flow rates have been established by using a contrast medium as stated on the catheter label at 37° centigrade for a 2 second injection. Since a number of variables (such as media viscosity and temperature) can produce different flow rates, this information should only be used as a guide.
- Reuse of single-use devices creates a potential risk of patient or user infections. Contamination of the device may lead to injury, illness or death of the patient.
- Reprocessing may compromise the integrity of the device and / or lead to device failure.


##### INSTRUCTIONS FOR USE

A. The following instructions outline the procedure for percutaneous entry utilizing the Seldinger technique. NOTE: Individual patient anatomy may preclude the utilization of this technique.


Puncture the vessel with the Seldinger needle.




1. Remove the cutting stylet; confirm flashback. Insert the distal end of guidewire through the needle cannula.

2. With the wire in position, remove needle cannula while maintaining compression at the puncture site.


3. Pass the dilator or sheath introducer assembly over guidewire into vessel. Remove dilator and insert catheter over guidewire.

4.

Remove guidewire, leaving catheter (or sheath introducer with catheter) in position.

5. B. Precautions to reduce clotting should be taken when any catheter is used. While in the vascular system, the catheter should always be filled with either sterile saline or contrast medium. Use of systemic heparinization should be considered.

C. Always confirm the position of the distal tip of the catheter using fluoroscopic guidance prior to injection using a pressure injector. Check for free backflow of blood, and inject a small amount of contrast medium by hand if necessary to confirm the catheter tip's position.

D. When using hydrophilic coated guidewires with ANGIODYNAMICS Angiographic Catheters, take special care to keep the guidewire soaking wet at all times, and follow the guidewire manufacturer's instructions for use. Consider using a smaller diameter guidewire when using 4 French catheters.

E. ANGIODYNAMICS Angiographic Catheters are designed for use with specific guidewire diameters. The recommended maximum guidewire diameter is specified on the catheter label. Guidewire use has been associated with a greater incidence of thrombus formation. Optimal guidewire size and judicious use are recommended.

F. ANGIODYNAMICS furnishes several of its catheters with a tip straightener. To straighten the catheter tip, hold the distal end of the tip straightener between the thumb and forefinger and pull the straightener up the catheter shaft until the tip is straight. Then place the catheter tip over the positioned guidewire, and introduce the catheter in the normal manner. To remove the tip straightener from the catheter, gently slide the straightener down the catheter shaft to the strain relief. Continue to gently pull until the straightener opens and slides off the catheter.

G. If resistance is encountered during catheter manipulation, determine the cause under fluoroscopy and take the necessary remedial action. If resistance is encountered when removing the guidewire from the catheter, simultaneously remove the catheter and the guidewire as a single unit under fluoroscopy to prevent potential vessel and product damage.

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