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IC 393 Rev. F
06-2014

**INDICATIONS:**

The AngioDynamics® VenaCure EVLT® NeverTouch® Procedure Kits are indicated for endovascular coagulation of the Great Saphenous Vein (GSV) in patients with superficial vein reflux, for the treatment of varicose veins and varicosities associated with superficial reflux of the Great Saphenous Vein (GSV), and for the treatment of incompetence and reflux of superficial veins of the lower extremity.

This product should be used only with lasers cleared for use in the treatment of varicose veins, varicosities with superficial reflux of the GSV, and in the treatment of incompetent refluxing veins in the superficial venous system in the lower limbs.

EXCLUSIONS:

- Patients with thrombus in the vein segment to be treated
- Patients with an aneurysmal section in the vein segment to be treated
- Patients with peripheral artery disease as determined by an Ankle-Brachial Index <0.9

POTENTIAL COMPLICATIONS:

The potential complications include but are not limited to the following: vessel perforation, thrombosis, pulmonary embolism, phlebitis, hematoma, infection, skin pigmentation alteration, neovascularization, paresthesia due to thermal damage of adjacent sensory nerves, anesthetic tumescence, non-target irradiation, hemorrhage, necrosis, skin burns and pain.

PRECAUTION:

- Prior to and during use, avoid damaging the fiber by striking, stressing, or excessive bending. Do not coil the fiber tighter than a diameter of 20cm. Clinical safety and effectiveness data is not available for other fiber tip designs and diameters.
- Prior to and during use, avoid bending the introducer sheath and dilator as this can cause kinks and damage.
- All persons in the treatment room must be wearing protective glasses with the proper rating for the wavelength being used.

WARNINGS:

- Treatment of a vein located close to the skin surface may result in skin burn.
- Paresthesia may occur from thermal damage to adjacent sensory nerves.
- Tissue not targeted for treatment must be protected from injury by direct and reflected laser energy. Use appropriate eye and protective wear for both patient and operating personnel.
- Reuse of single-use devices creates a potential risk of patient or user infections. Contamination of the device may lead to injury, illness, or death of the patient.
- Reprocessing may compromise the integrity of the device and/or lead to device failure.

INSTRUMENTS REQUIRED:

- Laser
- Safety eyewear

NEVERTOUCH KIT CONTENTS:

- NeverTouch Gold-Tipped Fiber with SMA 905 connector and Sheath-Lok® Fitting
- 4 French Trè-Sheath® Introducer and Dilator
- 0.018" and/or 0.035" guidewire
- Needle
- 5 French, 10cm sheath/dilator - (optional)

INDIVIDUALLY PACKAGED Trè-Sheath CONTENTS:

- 4 French Trè-Sheath Introducer and Dilator

PROCEDURE:

Carefully read all instructions and observe all Warnings prior to performing the procedure. Patient complications may occur if this is not done. This procedure can be performed in the physicians' office or as an outpatient treatment under local anesthesia and should be performed by a qualified physician who has received training in these techniques.

1. Carry out a physical examination.
2. Review Ultrasound vein diameter and depth measurements.
3. Mark the skin with a surgical pen to show the length of vein to be treated.
4. Set the laser up according to the provided "Setting up the Laser Procedure" operators manual.
5. Place the patient in reverse Trendelenburg position so as to engorge and distend the targeted treatment vein.
6. Using standard Seldinger technique gain access to the desired vein with the needle and corresponding guidewire.
7. If applicable, use the micro-access introducer to exchange for a larger diameter wire.
8. Advance the 4Fr Trè-Sheath introducer, over the wire, to the treatment location. If treating the Great Saphenous Vein, position sheath tip so that it is 1-2cm below the Sapheno-Femoral Junction. Verify sheath positioning using ultrasound guidance.
9. Remove the dilator and guidewire.
10. Flush the sheath with normal saline.
11. Properly attach the sterile, single use disposable NeverTouch laser fiber by removing the protective cap and completely screwing the fiber fitting to lasers' SMA 905 connector until tight and secure.
12. Insert the NeverTouch laser fiber into the sheath until the Fiber Stop Assembly comes in contact with the end of the Sheath Hub. (See below)

13. Remove Fiber Stop Assembly. Pull the sheath back and lock it to the Sheath-Lok Fitting.
14. Confirm location of the fiber using ultrasound guidance.
15. The NeverTouch fiber end should be 2 cm below the Sapheno-Femoral junction when treating the GSV.
16. Move the Sliding Sheath Gauge until flush with the insertion site, to mark the final position of the sheath.
17. Administer local tumescent anesthesia under ultrasound guidance as to ensure delivery of just enough anesthesia to achieve thermal protection.
18. Set the laser to continuous mode and adjust to desired power.
(810nm, 980nm Lasers: 10-14 Watts)
(1470nm Lasers: 5-7 Watts)
19. Put on laser safety glasses.
20. Place the laser in the Enable mode.
21. Check position of fiber end and sliding sheath gauge; adjust position if necessary.
22. Activate the laser by depressing the foot pedal while withdrawing the fiber and sheath together, at a rate that adequately delivers desired laser energy per centimeter.
(810nm, 980nm Lasers: 50-80 Joules per cm)
(1470nm Lasers: 30-50 Joules per cm)
Do not apply direct external hand pressure or force over the fiber tip during energy activation.
23. Cease operating the laser by removing foot from foot pedal when the fiber end is 2 – 3 cm from the access site as indicated by markers on the distal tip of the Trè-Sheath Introducer. The procedure is considered complete when the desired length of the targeted treatment vein is treated.

POST PROCEDURE:

Perform a second ultrasound to confirm that the vein is closed and that there is no more blood flow.

WARRANTY

AngioDynamics, Inc. warrants that reasonable care has been used in the design and manufacture of this instrument. This warranty is in lieu of and excludes all other representations and warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose. Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond AngioDynamics, Inc.'s control directly affect the instrument and the results obtained from its use. AngioDynamics, Inc.'s obligation under this warranty is limited to the repair or replacement of this instrument and AngioDynamics, Inc. shall not be liable for any incidental or consequential losses, damages or expenses directly or indirectly arising from the use of this instrument. AngioDynamics, Inc. neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. AngioDynamics, Inc. assumes no liability with respect to instruments reused, reprocessed or resterilized, modified or altered in any way, and makes no warranties or representations, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.

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Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

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*AngioDynamics, the AngioDynamics logo, NeverTouch, VenaCure EVLT, Trè-Sheath, and Sheath-Lok are trademarks and/or registered trademarks of AngioDynamics Inc., an affiliate or a subsidiary. May be covered by or for use under one or more of the following U.S. Patents: 5,714,360; 6,986,766; 7,033,347; 7,559,329; 7,483,457; 7,907,643; 8,413,664; and other U.S. and Foreign counterparts.
© AngioDynamics, Inc. 2014 All rights reserved

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© AngioDynamics, Inc. 2014, minden jog fenntartva.
© AngioDynamics, Inc. 2014. Všechna práva vyhrazena.
© AngioDynamics, Inc. 2014. Wszystkie prawa zastrzeżone.
© AngioDynamics, Inc. 2014 Med enerett.
©AngioDynamics, Inc. 2014 Tüm hakları saklıdır.
©AngioDynamics, Inc., 2014. Visos teisės saugomos.