LIFEGUARD™

s a f e t y n e e d l e

INSTRUCTIONS FOR USE

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MODE D'EMPLOI

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ANGIODYNAMICS®

ANGIODYNAMICS®

TABLE OF CONTENTS

en	INSTRUCTIONS FOR USE	1
du	GEBRUIKSAANWIJZINGEN	9
es	INSTRUCCIONES DE USO	18
it	ISTRUZIONI PER L'USO	27
fr	MODE D'EMPLOI	36
gm	GEBRAUCHSANWEISUNG	45

PRODUCT DESCRIPTION:

The LifeGuard™ Safety Infusion Set is lipid compatible, latex free and is constructed with Non-DEHP PVC.

The LifeGuard™ Safety Infusion Set is designed to reduce needlestick injuries when used according to the enclosed directions. Carefully read and follow all instructions prior to use of the device.

The LifeGuard™ Safety Infusion Set is a 90° non-coring needle composed of wings, tubing, a universal luer connector, optional Needleless Y-site, tubing clamp(s), needlestick guard and needle trap that encapsulates the needle when de-accessing the port.

The optional Needleless Y-Site has been incorporated into the LifeGuard™ Safety Infusion Sets for your convenience. Please refer to the manufacturer's directions for use, which is included in the packaging of those needle product codes.

Inner package contents sterile and non-pyrogenic, if not opened or damaged. Sterilized by Ethylene Oxide (EtO).

STERILE EO

INDICATIONS FOR USE:

The LifeGuard™ Safety Infusion Set is used to access implanted vascular ports to administer fluids and/or to withdraw blood.

The LifeGuard™ Safety Infusion Set facilitates the safe removal of the needle by encapsulating the needle during vascular port de-accessing to help prevent needlestick injuries.

The LifeGuard™ 19 and 20 Gauge non Y-Site needles may be used with the AngioDynamics® Smart Port™ CT implantable ports for power injection of contrast media. Maximum pressure recommended for power injection of contrast media is 300 psi and the maximum flow rate is 5 ml/sec. For information and

instructions regarding the power injection procedure, please refer to the Instructions for Use Manual for the Smart Port™ CT Power Injectable Implantable Port Systems.

MRI INFORMATION:

The LifeGuard Safety Infusion Set was determined to be MR-conditional according to the terminology specified in the American Society for Testing and Materials (ASTM) International, Designation: F2503-05. Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment. ASTM International, 100 Barr Harbor Drive, PO Box C700, West Conshohocken, Pennsylvania, 2005.

Non-clinical testing demonstrated that the LifeGuard Safety Infusion Set is MR Conditional. A patient with this device can be scanned safely immediately after placement under the following conditions:

- Static magnetic field of 3-Tesla or less
- Spatial gradient magnetic field of 720-Gauss/cm or less

MRI-Related Heating

In non-clinical testing, the LifeGuard Safety Infusion Set produced the following temperature rise during MRI performed for 15-min in the 3-Tesla (3-Tesla/128-MHz, Excite, Software G3.0-052B, General Electric Healthcare, Milwaukee, WI) MR system:

Highest temperature change +1.7° C

Therefore, the MRI-related heating experiments for the LifeGuard Safety Infusion Set at 3-Tesla using a transmit/receive RF body coil at an MR system reported whole body averaged

SAR of 3.0-W/kg (i.e., associated with a calorimetry measured value of 2.8-W/kg) indicated that the greatest amount of heating that occurred in association with these specific conditions was equal to or less than +1.7° C.

Artifact Information

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the LifeGuard Safety Infusion Set. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary.

Pulse Sequence	T1-SE	T1-SE	GRE	GRE
Signal Void Size	1,299-mm ²	1,130-mm ²	4,987-mm ²	4,081-mm ²
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular

CONTRAINDICATIONS:

Do not use with patients showing signs or symptoms of infection at the insertion site.

Do not use with patients who have exhibited prior intolerance to the materials of construction.

WARNINGS:

Review all fluids to be used with this product for any incompatibility with plastics and plastic adhesives.

A 10 ml syringe or larger is recommended.

Maximum pressure recommended for power injection of contrast media with the Smart Port™ CT implantable ports is 300 psi and a maximum flow rate of 5ml/sec. Only LifeGuard™ 19 and 20 gauge non Y-Site needle models may be utilized for power injection with Smart Port™ CT.

Do not use the device if the package has been opened or damaged.

Clamp infusion set when not in use.

Store in a cool dry place.

CAUTIONS:

Federal (USA) law restricts this device to sale by or on the order of a physician.

RX only

SINGLE USE ONLY: Non-Pyrogenic, Non-toxic. Contents sterile unless package is opened or damaged. Do not resterilize.

Carefully read and follow all instructions prior to use of the device.

Follow CDC guidelines, institutional guidelines or INS and ONS guidelines for infusion set use.

Always maintain universal precautions and utilize aseptic technique throughout insertion care and maintenance procedures.

Application of excessive force to needle while accessing port may result in damage to needle and port.

LEAKAGE OF FLUID MAY OCCUR WHEN DISCONNECTING COMPONENTS.

KEEP HAND / FINGERS PROTECTED FROM NEEDLE TIP AT ALL TIMES BY FOLLOWING THE INSTRUCTIONS FOR USE.

DO NOT ATTEMPT TO OVERRIDE THE SAFETY FEATURE.

INSTRUCTIONS:

PROCEDURE FOR PORT ACCESS:

Note: Aseptic technique should be maintained during port access.

Note: Follow Institutional Universal Precautions.

1. Gather supplies necessary for port access.
2. Place mask on face.
3. Open and drop all supplies onto sterile field.
4. Don sterile gloves.
5. Cleanse area around port, following your institutional protocol.
6. Remove luer cap(s) and attach injection cap(s) used by your facility.

Note: Prior to infusing, ensure that all injection caps are securely connected.

7. Prime infusion set following your institutional protocol.

INSERTION:

1. Grasp the flexible wings with thumb and middle finger, placing your index finger on top of the needlestick guard (see Figure A-page 8).
2. Remove protective sleeve on needle.
3. Insert needle at a 90° angle to the port.
4. Advance needle through the skin and the port septum until it contacts the bottom of the port reservoir.

Note: Needle placement should be confirmed by aspiration of 3-5 ml of blood.

Dual port access: When accessing septum #1 of the port, place the first needle with the needlestick guard facing laterally to the patient. When accessing septum #2 of the port, place the second needle with the needlestick guard facing medially to the patient.

5. Secure the device per institutional protocol.
6. Using turbulent “push-pause” technique, flush port per institutional protocol maintaining positive pressure.

Note: A 10 ml syringe or larger is recommended.

REMOVAL:

7. Using turbulent “push-pause” technique, flush the non-coring needle per institutional protocol, maintaining positive pressure.
8. Remove dressing.

Tip:

Grasp the upper right-hand corner of the dressing and the lower left-hand corner of the dressing and pull and stretch outward and upward. Care should be used to ensure the needle does not become dislodged from the patient’s port. While holding the needle securely against the skin, continue removing the dressing until the wings are exposed.

Note: Alcohol may be used to assist with dressing removal. You may continue to deaccess the needle with the dressing still covering the needle trap and needlestick guard.

9. Raise the needle trap to a 90° angle (see Figure B-page 8).
10. Using your non-dominant hand, grasp the needlestick guard with your thumb and index finger. Gently push the needlestick guard down against the patient's port.
11. While securely holding the needlestick guard in place, use your dominant hand to grasp the flexible wings and pull upward until the needle is completely encapsulated in the needle trap (see Figure C-page 8).

Note: The needle trap allows for visual confirmation that the needle is fully encapsulated and safe. Additionally, you will hear it lock into the safe position.

12. Flip the needlestick guard towards the needle trap (See Figure D-page 8).
13. Properly dispose in Sharps Container.



INJECTION CAP CHANGE RECOMMENDATION:

1. Change injection caps per institutional guidelines or when cap has been removed for any reason, appears damaged, is leaking, or if there is blood residue on the cap.

INJECTION CAP CHANGE PROCEDURE:

1. Wash hands thoroughly with soap and warm water.
2. Don clean gloves.
3. Ensure tubing clamp is closed prior to disconnecting the injection cap.
4. Open the sterile injection cap package.
5. Purge air from device by flushing with appropriate solution.
6. Clean the tubing's luer lock hub on the outside with povidone-iodine or swab with alcohol. Allow to dry.
7. Remove the protective piece from the new injection cap, being careful not to touch the inner portion of the new

- injection cap or infusion set luer.
8. Attach the new injection cap to the luer in a clockwise manner. DO NOT OVER TIGHTEN.

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
← FLEXIBLE FLÜGEL

DEVICE COLOR CODING

Brown = 19 Gauge Needle
Yellow = 20 Gauge Needle Black = 22 Gauge Needle

Keep Dry

Keep away from sunlight.

Do not use if package is damaged.
Niet gebruiken als de verpakking is beschadigd.
No utilizar si el envase está danado.
Non utilizzare se la confezione è danneggiata.
Ne pas utiliser le dispositif si l'emballage est endommagé.
Nicht verwenden, wenn die Verpackung beschädigt ist.

LifeGuard is a trademark of AngioDynamics, Inc.

US Patent #6,676,633

NON-PYROGENIC

STERILE **EO** STERILIZED WITH
ETHYLENE OXIDE

REFER TO ENCLOSED
INSTRUCTIONS FOR
PROPER USE OF THIS DEVICE

SINGLE USE ONLY. DO NOT
RESTERILIZE. STERILE IF
PACKAGE IS NOT OPENED OR
DAMAGED

CAUTION: Federal (U.S.A.) Law restricts this device to sale by or on the order of a physician.

ANGIODYNAMICS®

Manufactured by:
AngioDynamics, Inc.
One Horizon Way
Manchester, GA 31816
P: 518-798-1215 or
800-772-6446

www.angiodynamics.com