
RITA® StarBurst® Access System, Soft Tissue


INDICATIONS:
To provide access for delivery of the RITA® StarBurst® devices, fluids and material such as saline to interstitial tissue.

INSTRUCTIONS FOR USE:


Placement of the Introducer
• When using the introducer to place the StarBurst® device, the introducer is used to place the StarBurst® device in the target area with the stylet tip positioned where the tip of the StarBurst® device will be.

Placement of the Device Inside the Introducer
• Ensure that the introducer is not bent.
• Place the StarBurst® device **slowly** into the introducer.
• The device should go through the introducer **very easily** if the introducer is straight. If resistance is felt during insertion of the device, stop and make sure the introducer is straight before proceeding. If the device is forced through a bent introducer it could cut into the introducer.

• The unutilized portion of the trocar of the StarBurst® device should be **subdistal** to the introducer until the hub of the StarBurst® device touches the proximal end of the introducer.
• If inserting a StarBurst® device only partially into an introducer, insert the StarBurst® device **upward** until reaching the appropriate marking on the StarBurst® trocar according to the table below.

Introducer/StarBurst® Device Compatibility


StarBurst® Devices	10 cm (700-102636)	13 cm (700-102637)	23 cm (700-102638)
XL (10 cm)	Not compatible	Not compatible	Not compatible
XL Semi-Flex (25 cm)	To 12 cm mark	To 15 cm mark	To Hub
XL (15 cm)	To 12 cm mark	To Hub	Not compatible
XL (25 cm)	To 12 cm mark	To 15 cm mark	To Hub
XLie (12 cm)	To Hub	Not compatible	Not compatible
XLie (25 cm)	To 12 cm mark	To 15 cm mark	To Hub
XLie SemiFlex	Not compatible	Not compatible	Not compatible
Talon (15 cm)	To 14 cm mark	To 17 cm mark	Not compatible
Talon Semi-Flex (25 cm)	Not compatible	Not compatible	Not compatible

Note: If there is any resistance on deployment, ensure that the StarBurst® device is properly inserted into introducer to the point indicated above.

Precautions
• **CAUTION:** Do not attach anything abrasive (e.g. metal clamps, etc.) to the introducer and/or Device. This may damage the introducer, which could contribute to patient injury.
• **CAUTION:** Care must be taken when inserting the introducer into hard and/or abrasive tissue. This may damage the introducer, which could contribute to patient injury.

• **CAUTION:** Care must be taken when inserting the introducer into hard and/or abrasive tissue. This may damage the introducer, which could contribute to patient injury.

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STERILE/LE Method of Sterilization: Ethylene Oxide

CE 008 CE Mark: Indicates conformance with European Council Directive 93/42/EEC

AngioDynamics, Inc. One Horizon Way, Manchester, GA 31816
AngioDynamics UK Ltd. Building 2000, Beach Drive, IO Cambridge, Waterbeach, Cambridgeshire, CB25 9TE, United Kingdom
Tel: +1-800-772-6448, +1-518-798-1215, +1-518-798-1360
Fax: +44 (0)122-372-9361, +44 (0)122-372-9329

RITA and StarBurst are trademarks of AngioDynamics, Inc. For warranty information, refer to the RF Generator User's Guide.

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