VenaCureEVLT™ Procedure Kit With Spotlight OPS™ Sheath For Endovenous Laser Treatment

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Manufactured by:

AngioDynamics UK Limited, Building 2000, Beach Drive, IQ Cambridge, Waterbeach,
Cambridge, CB25 9TE, United Kingdom

Reorder from:

US/Canada: Tel: +1 800 772 6446 Fax: +1 518 798 1360

Rest of World: Contact your local AngioDynamics Distributor or Tel: +44 1223 729300

www.venacureevlt.com or www.angiodynamics.com

Covered by or for use under one or more of the following US Patents:
6,398,777; 6,769,433; 6,752,803; 6,258,084; 6,981,971; 7,396,355,
their foreign counterparts, and other patents pending

Rx Only (USA)

PP/1279 Iss 3

NOTE: Read this manual and the Laser Operator's manual thoroughly prior to using the VenaCureEVL™ procedure kit. Observe all warnings, precautions and cautions noted throughout these documents. Failure to do so may result in patient complications.

ATTENTION: **SINGLE USE ONLY**

ATTENTION: Insure expiration dates are still valid.

ATTENTION: Inspect the sealed packages before opening. If seals are broken or the packages are damaged, treat as non-sterile and discard.

Disposal: Used components should be disposed of in accordance with local regulations or institutional policy.

CAUTION: Intended for use only by fully trained physicians.
Federal Law (USA) restricts this device to sale by or on the order of a physician.

VenaCureEVL™ Procedure kit contents

- PACK A
 - Percutaneous entry needle
 - Double ended Guidewire
- PACK B
 - Optical fiber with distance markers (Site Marks™)
 - Introducer sheath with fiber clamp and entry point marker.

Additional equipment required but not supplied

- Diode Laser with suitable Laser Safety glasses
- Duplex Ultrasound

In addition to the above a variety of consumable and disposable items will be required such as local anesthetic, drapes, syringes and needles. These will vary according to physician preference.

Indications

The Laser and VenaCureEVL™ Procedure Kit are intended for use in the treatment of varicose veins and varicosities associated with superficial vein reflux of the Great Saphenous vein, and with veins in the lower limbs with superficial reflux.

The D15 Plus, D30 Plus, DELTA-15 and DELTA-30 lasers and VenaCureEVL™ Kits are indicated for treatment of incompetent refluxing veins in the superficial venous system.

Contraindications

- Patients with thrombus in the vein segment to be treated.
- Patients with an aneurysmal section in the vein segment to be treated.
- Patients with peripheral arterial disease as determined by an Ankle-Brachial Index < 0.9.
- Patients with an inability to ambulate.
- Patients with deep vein thrombosis (DVT).
- Patients who are pregnant or breast feeding.
- Patients in general poor health.
- Other contraindications may be raised by the individual physician at the time of treatment.

Extremely tortuous vein segments may require treatment by alternative techniques (phlebectomy, sclerotherapy).

Diomed Equipment Handling Requirements

Diode Laser

- Prepare for use in accordance with the Laser Operator Manual.

VenaCureEVL™ Procedure Kit (PACKS A & B)

- Using sterile technique, remove both Packs (A & B) from the outer packing.
- Packs A & B may be placed in the sterile field if desired.

VenaCureEVL™ Procedure Kit (PACK A)

- Using sterile technique open the pack, remove all contents, and inspect for damage. DO NOT USE if any component is damaged.
- If damage is seen, contact Customer Service or your local distributor.

VenaCureEVL™ Procedure Kit (PACK B)

- Do NOT open until access has been made to the vein using the contents of PACK A
 - Using sterile technique open the pack, remove all contents, and inspect for damage. DO NOT USE if any component is damaged.
- If damage is seen, contact Customer Service or your local distributor.

VenaCureEVL™ Laser Fiber

- When required, remove the sterile fiber from the package, keeping it coiled.
- Remove the securing ties and gently uncoil the fiber. Do not remove the SMA connector cover at this time.

Instructions for Use

1. Prepare and drape the leg using standard practice.
2. Use sterile technique to remove and open Pack A.
3. Place the patient in the reverse Trendelenburg position and use the percutaneous needle to gain access to the vein (under local anesthesia and ultrasound guidance).

NOTE: USE ULTRASOUND GUIDANCE THROUGHOUT THE PROCEDURE to ensure that devices are correctly located.

4. Advance the guidewire through the needle until it reaches the junction with the deep vein system.
5. Remove and discard the percutaneous needle.
6. Use sterile technique to remove and open Pack B.
7. Insert the VenaCureEVL™ Spotlight OPS sheath over the guidewire until it reaches the junction with the deep vein system.

CAUTION: DO NOT USE EXCESSIVE FORCE to introduce the sheath. Excessive force may cause damage to, or separation of the sheath.

8. Remove and discard the guidewire.
9. Advance the VenaCureEVL™ optical fiber through the Spotlight OPS sheath until the distal SiteMark on the optical fiber is coincident with the compression clamp (back of the sheath).

CAUTION: DO NOT USE EXCESSIVE FORCE to introduce the optical fiber. Use ultrasound to confirm that the Spotlight OPS sheath is correctly positioned if strong resistance is encountered.

10. Use ultrasound guidance in conjunction with the echogenic tip on the Spotlight OPS sheath to confirm that the tip of the sheath is 1-2cm below the junction with the deep vein system.
11. While keeping the optical fiber in its location, withdraw the sheath until the compression clamp is coincident with the proximal SiteMark on the optical fiber.
12. Tighten the compression clamp.
13. Slide the Entry Point Marker (EPM) to the tissue access site to record the length of vein treated.
14. Place the patient in the Trendelenburg position and administer tumescent anesthesia along the entire length of the vein being treated.
15. Attach the optical fiber to the laser.

WARNING: LASER PROTECTIVE EYEWARE MUST BE WORN BY EVERYONE IN THE TREATMENT ROOM.

16. Use trans-illumination of the red aiming beam through the skin, and ultrasound to make final confirmation that the tip of the optical fiber is more than 1cm from the junction with the deep vein system.

17. Activate the laser and using the distance markers along the length of the Spotlight OPS sheath, withdraw the sheath (with optical fiber) at a rate consistent with the venous anatomy.

CAUTION: STOP THE LASER FIRING BEFORE THE OPTICAL FIBER IS WITHDRAWN THROUGH THE ACCESS SITE.

18. Instruct the patient to ambulate.

WARNING: INSTRUCT THE PATIENT TO AVOID HOT BATHS AND VIGOROUS ACTIVITY FOR 7 DAYS FOLLOWING VenaCureEVL™.

19. Dress the patient's leg and schedule subsequent follow-up according to physician preference.

General Guidance

Potential complications include, but are not limited to the following: vessel perforation, thrombosis, pulmonary embolism, phlebitis, hematoma, infection, paresthesia, skin burns, and thrombophlebitis.

WARNING: TREATMENT OF A VEIN LOCATED CLOSE TO THE SKIN SURFACE MAY RESULT IN A SKIN BURN.

WARNING: PARESTHESIA MAY OCCUR FROM THERMAL DAMAGE TO ADJACENT SENSORY NERVES.

WARNING: TISSUE NOT TARGETED FOR TREATMENT MUST BE PROTECTED FROM INJURY BY DIRECT AND REFLECTED LASER ENERGY WITH APPROPRIATE EYE AND PROTECTIVE WEAR FOR ANY PERSON PRESENT IN THE OPERATING ROOM.

PRECAUTION: PRIOR TO AND DURING USE, AVOID DAMAGING THE FIBER BY STRIKING, STRESSING OR EXCESSIVE BENDING. DO NOT COIL THE FIBER TIGHTER THAN A RADIUS OF 60mm.

PRECAUTION: THE POSITIONS OF THE SITEMARKS ON THE VenaCureEVL™ FIBER HAVE BEEN MATCHED TO THE INTRODUCER SHEATH PROVIDED IN THE VenaCureEVL™ PROCEDURE KIT. ALTERNATIVE SHEATHS MUST NOT BE SUBSTITUTED.

PRECAUTION: PRIOR TO AND DURING USE, AVOID TIGHT BENDING THE INTRODUCER SHEATH AS THIS CAN CAUSE KINKS AND DAMAGE. DO NOT COIL THE SHEATH TIGHTER THAN A RADIUS OF 60mm.

PRECAUTION: THE ENTRY NEEDLE PROVIDED WITH THIS KIT IS MATCHED TO THE SIZE OF THE GUIDEWIRE AND SHOULD NOT BE SUBSTITUTED.

PRECAUTION: DO NOT TIGHTEN THE COMPRESSION CLAMP ON SHEATH UNTIL FIBER IS IN POSITION.

NOTE: PATIENTS ON ANTICOAGULANT THERAPY ARE STILL ELIGIBLE FOR TREATMENT. MOST PHYSICIANS WILL ONLY TREAT PATIENTS WITH AN I.N.R. (INTERNATIONAL NORMALISED RATIO) AT THE LOWER THERAPEUTIC LIMITS.

NOTE: IT IS SOMETIMES HELPFUL TO NICK THE SKIN WITH A SCALPEL AT THE PUNCTURE SITE EITHER TO AID INSERTION OF THE NEEDLE, OR AFTER GUIDEWIRE INSERTION.

VenaCureEVL™ can be performed under local anesthesia in the physician's office or as an outpatient treatment.

VenaCureEVL™ must only be performed by a qualified physician who has received specialized training in the techniques described.

Storage

- The VenaCureEVL™ Procedure Kit should be stored in a cool, dry and dark area.
- Avoid extended exposure to direct sunlight and extremes of temperature, since this may cause degradation or discoloration of the packaging and contents.

Do not open the package until ready to use the VenaCureEVL™ Procedure Kit.

Symbols used on packaging

Single Use Only

Refer to Operator Instructions

Batch Number

Sterile by Ethylene Oxide

Expiration Date

Product Re-order code

Do Not Use If Packaging Is Damaged

No Latex

Usable Sheath Length

This Way Up

Store Out of Direct Sunlight

Keep Dry

Fragile

Note: VenaCureEVL™, Spotlight OPS and SiteMarks are a registered trademark of Angiodynamics.

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