PATIENT EDUCATION

## **IMPLANTABLE PORT SYSTEMS**

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INFORMATION DU PATIENT

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**RITA Medical Systems**

# RITA Medical Systems

## IMPLANTABLE PORT SYSTEMS

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# PATIENT EDUCATION

*Dear Patient,*

*Your doctor has recommended that you receive an implantable port. To introduce you to the implantable port, we have prepared this brief pamphlet. As with any medical treatment, there are risks associated with the use of an implantable port, but this pamphlet can provide only general information. We therefore urge you to discuss any questions you have about the device or your treatment with your doctor or nurse. You should follow your doctor's advice and instructions concerning your implantable port and medical treatment. For your convenience, we have included a patient identification card.*

*The Staff at RITA*

THE IMPLANTABLE PORT ADVANTAGE	<b>2</b>
DESCRIPTION OF THE IMPLANTABLE PORT	<b>2</b>
PLACEMENT OF THE IMPLANTABLE PORT	<b>3</b>
USE OF THE IMPLANTABLE PORT	<b>3</b>
POTENTIAL PROBLEMS WITH THE IMPLANTABLE PORT	<b>4</b>
CARE OF THE IMPLANTABLE PORT	<b>6</b>
QUESTIONS REGARDING YOUR IMPLANTABLE PORT	<b>6</b>
PATIENT IDENTIFICATION CARD	<b>8</b>


## The Implantable Port Advantage

RITA Medical Systems, Inc. (RITA) manufactures implantable access devices designed to help people continue their normal activities while receiving medical treatments.

The unique feature of an implantable port is that it is placed completely under the skin. This placement allows almost complete freedom of movement for the patient. When the port is in place, repeated access to the bloodstream or a specific body site to administer therapy or withdraw blood can be accomplished with greater ease and less discomfort. In addition, the port requires minimal care.

## Description of the Implantable Port

The port is a small plastic or metal disc, about 2-3 cm in diameter (the size of a half dollar), with a slightly raised rubber injection site called the port septum. Since the septum is raised, it is easily identified from the surface of your skin, making injections quick and easy. The septum is made of a silicone material which reseals itself after each use. The material allows the septum to be punctured many hundreds of times with a special needle and yet retain its integrity.

Attached to the base of the port is a narrow, flexible tube, called a catheter. The catheter is usually inserted into a large blood vessel and delivers your therapy from the port body into your bloodstream.


Implantable ports are available in a variety of sizes and models. For example, dual ports have two separate septa contained

within one port body, which allows the port to be used for two different purposes simultaneously. Your doctor will decide which model is appropriate for you, based on your individual needs.

## Placement of the Implantable Port

Although each person's situation is unique, the port is usually inserted during a brief surgical procedure performed under local anesthesia. You may even have your port inserted in a Day Surgery or Ambulatory Surgery Unit and go home shortly afterwards.

During the procedure, your surgeon will make two small incisions. The catheter is placed into the selected blood vessel or body site through one incision and then tunneled under the skin to the other incision site. Here, a “pocket” is formed under the skin to

hold the port. The catheter is attached to the port and the port is then sutured in place securely within the pocket. An x-ray will be performed to verify that the catheter is in the proper position in the blood vessel. A few stitches are required to close the pocket once the port is in place.

Bandages may be taped over the incisions for the first few post-operative days. The skin overlying the port will be swollen and tender, but this irritation should go away as the incisions heal. The port may be used immediately if necessary, or your doctor may wait until the swelling has decreased. During these first few days after surgery, it is important that you avoid any heavy exertion or strenuous activities.

## Use of the Implantable Port

The most frequent use of a port is for venous access. “Venous” refers to certain blood vessels in the body. “Access” in this context, means a way of entering these vessels.

Using a special anti-coring needle, your doctor or nurse can go through the skin overlying the port and into the center of the septum with one quick needle stick. These needles can then be used for a variety of purposes, such as administration of intravenous fluids, blood products, and medications, as well as withdrawal of blood samples for laboratory tests.

Injections of medications can be done over a very short period of time (called a “bolus”) or over a prolonged period of time (called a “continuous infusion”). Continuous infusions may last anywhere from one hour up to several days or longer and may require the use of a portable infusion pump to control the rate of the infusion. If your medication is to be given over a long period of time, the needle will stay in place, covered with a sterile dressing, and be changed at intervals as decided by your doctor. Once the needle is removed, the dressing is no longer needed.

## Potential Problems with the Implantable Port

The use and maintenance of long term vascular access devices may cause possible problems or complications. Listed below are the signals associated with the more common complications. If any of these symptoms appear, or any other unusual symptoms, it is important to contact your doctor or nurse immediately.

<u>What you are experiencing?</u>	<u>Possible Reason</u>
Inability to withdraw blood, or inability to flush the port and catheter using normal pressure.	Catheter Occlusion or blockage
Aching discomfort, swelling or pain in the shoulder, neck or arm.	Vein Thrombosis
Pain or redness over or around the port site, fever, or any unusual drainage from the incision site.	Infection of the port pocket or at the needle insertion site
Burning sensation upon infusion of medication	Extravasation/ Infiltration—the leakage of drug into the tissue surrounding the port

Other complications include, but are not limited to:

- Catheter Pinch-Off- the catheter being “pinched” between the clavicle and first rib-can lead to Catheter Shear, which results in breakage of the catheter
- Device Rotation or Dislodgement can occur if the port is not well sutured in place or with excessive manipulation
- Catheter Migration can occur spontaneously and result in the catheter being positioned outside of the intended vessel
- Catheter Disconnection is the separation of the catheter from the port
- Skin Erosion can occur over the port site
- Complications associated with surgery
- Rejection of the device by the patient’s body

## Care of the Implantable Port

Because the port is implanted completely beneath the skin, it requires minimal care. If permitted by your doctor and subject to his or her instructions, you may bathe, swim and exercise as you always have, without worry about damaging the port. Once the skin has healed after surgery, there is no need to cover your port site with a dressing when the port is not accessed with a needle. In between treatments, you may not even remember that you have an implantable port, since it will be visible only as a small bump on your skin.

Periodically, your port will need to be “flushed” to ensure that it remains in working order. Your doctor or nurse will let you know how often this needs to be done and who will perform this flushing procedure. In addition, you should examine your port site regularly. Once the initial expected tenderness and swelling during the post-operative period has subsided, you should notify your doctor or nurse if you notice any additional swelling, discoloration (bruising), redness, tenderness, pain, or drainage at the incision sites, the needle insertion site, or the port pocket, or if you begin to run a fever or experience any other unusual symptoms.

## QUESTIONS REGARDING YOUR IMPLANTABLE PORT

### Why do I need an implantable port?

Your doctor has recommended the implantable port for you because it is a reliable route for the administration of therapy for your specific medical needs. The port allows easy repeated access to the bloodstream or a specific body site, which in turn causes less vessel or site damage.

## **Will I have to be stuck with a needle and will it be painful?**

You will need to have a needle inserted into the port septum to administer therapy or withdraw blood, but the access site will be easy to locate. You will notice some swelling and tenderness at the implantation site immediately following the port insertion. This may cause some pain and discomfort when your nurse or doctor inserts a needle into the port. However, over time, the initial tenderness will subside and needle insertions should become less painful.

## **Will my normal activities be affected by my port?**

In the immediate period following the port placement procedure, you should avoid strenuous activities. Once your incisions heal, your normal activities can be resumed. With your doctor's approval, you may bathe, swim and exercise while your port is not accessed with a needle.

## **Will my insurance company pay for the port?**

This varies from insurance company to insurance company. In general, most companies provide coverage for the placement and care of a port.

## **Will my port be affected by x-rays, MRI, or other diagnostic tests my doctor may order?**

There are no components in an implantable port that could be affected by various methods of imaging. Your physician or nurse may call our Customer Service Department if any questions arise.

## **How long can the port remain implanted and can it be removed once it is no longer needed?**

The port may remain in place for extensive periods as long as your physician feels your medical treatment requires it. Once your physician has determined that your port is no longer needed, it can be removed during a brief surgical procedure.

## **Will my port be visible?**

The port is implanted completely under the skin. During the immediate period following implantation, there may be some redness and swelling. This irritation should disappear in a few weeks. There will be a small raised area over the port which can be felt and may be visible. There also will be a small scar located near the port, which should also diminish in visibility over time.

## **Patient Identification Card**

It is very important for people who have implanted devices to always carry an identification card with pertinent medical information. This is particularly important in the event of an emergency, so that the medical personnel will be aware that you have a port. For your convenience, we have included such a card.

Prior to leaving the hospital, please have the surgeon who inserted your port, or one of the nursing staff, help you fill in the required information. Remember to have this card with you when you visit your doctor, hospital, clinic, or dentist so that the information on your port will be immediately available.

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CAUTION: Federal (U.S.A.) Law restricts this