



Radiofrequency Ablation Electrodes

StarBurst* Talon

INFUSION RFA ELECTRODES

BENEFITS

- Integrated disposable main cable and tubing allows for easy set-up
- Solid tip allows easier penetration of tissue
- Creates spherical ablations up to 4 cm
- Dynamic, real-time temperature readouts
- Infusion-based system increases conductivity of ablation zone
- Needle track ablation
- Used for percutaneous, laparoscopic and open cases


SPECIFICATIONS

- Four active arrays plus active trocar tip
- Thermocouple in tip of each active array


StarBurst Talon

INFUSION RFA ELECTRODES

4 cm Ablation - 9 minutes^{††}



^{††}Hold time at target temperature (set timer)



Product Name	Length	Part #	Outer Diameter (O.D.)
StarBurst Talon Electrode	15 cm	700-102847	14 gauge/6.4 French
StarBurst Talon Electrode	25 cm	700-102846	14 gauge/6.4 French
StarBurst Talon Semi-Flex Electrode	25 cm	700-102845	14 gauge/6.4 French (rigid portion)

Model 1500X RF Generator & IntelliFlow Pump

The Model 1500X RF Generator is designed specifically for use with RITA[®] electro-surgical devices. It is the latest radiofrequency generation system that features technological advances including software upgrade capabilities, potential for 250 watts of power, and three flexible serial ports. The IntelliFlow Pump is used in conjunction with the 1500X RF Generator for the ablation of soft tissue and offers the ultimate in ease of use, while delivering a precise amount of infused saline.



IMPORTANT RISK INFORMATION

INDICATIONS FOR USE: To be used in conjunction with the RITA RF Generator and IntelliFlow Infusion pump for the ablation of soft tissue.

WARNINGS: Do not attach anything (i.e., clamps, etc.) to the Device. This may damage the insulation, which could contribute to patient injury. Do not twist or exert high forces on the device while it is deployed in the tissue. This may cause the needles to break and remain in the tissue. Do not remove the Device without ensuring that the needles are fully retracted within the trocar. This may cause the needles to break and remain in the tissue. If the device is being used in a laparoscopic procedure, care must be taken to avoid a gas embolism. For SemiFlex Talons: over bending of the trocar to a radius smaller than 5 cm beyond a 90° curvature and/or kinking the device can damage the trocar and cause

patient injury. Patients with peripheral vascular deficiency are at increased risk of thermal injury from Dispersive electrodes. Patients with frail skin are at increased risk of skin damage from the adhesive on the Dispersive pads. Do not use metal introducers that do not have insulation. RF energy can be transmitted from the electrode through the un-insulated metal introducers to the patient, causing inadvertent burns. The Talon Electro-surgical Device is MR UNSAFE.

PRECAUTIONS: Having RF power on at the same time as infusion using a method different from these instructions may alter the path of the electrical energy away from target tissues. Do not bend or kink the trocar or the needles. This may cause damage and result in a non-functional device. Keep the voltage/power as low as possible to achieve the desired end effect. If the

device is being used in a laparoscopic procedure, activation of the device when not in contact with target tissue may cause capacitive coupling. Reuse of single-use devices creates a potential risk of patient or user infections. Contamination of the device may lead to injury, illness or death of the patient. Reprocessing may compromise the integrity of the device and/or lead to device failure.

Refer to individual product IFUs and/or User Manual to see full Warnings, Precautions, Possible Adverse Effects and Contraindications. Observe all instructions prior to use. Failure to do so may result in patient complications.

CAUTION: Federal (USA) law restricts these devices to sales by or on the order of a physician.



USA > 14 Plaza Drive, Latham, NY 12110 > tel: 800-772-6446 or 518-798-1215 > fax: 518-798-1360
International > Haaksbergweg 75 (Margrietoren), 1101 BR, Amsterdam Z-O > The Netherlands
tel: +31 (0)20 753 2949 > fax: +31 (0)20 753 2939

www.angiodynamics.com

*AngioDynamics, the AngioDynamics logo, StarBurst, the StarBurst logo and RITA are trademarks and/or registered trademarks of AngioDynamics, Inc., an affiliate or a subsidiary.

© 2013 AngioDynamics, Inc. MLC 383 US Rev A