**RITA® StarBurst® XLi-enhanced and  
StarBurst® XLie SemiFlex  
Electrosurgical Device**

**INSTRUCTIONS FOR USE**

**BRUGSANVISNING**

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**MODE D'EMPLOI**

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**ANGIODYNAMICS®**

# ANGIODYNAMICS®

## TABLE OF CONTENTS

<b>[en]</b>	<b>INSTRUCTIONS FOR USE</b>	<b>1</b>
<b>[da]</b>	<b>BRUGSANVISNING</b>	<b>8</b>
<b>[de]</b>	<b>GEBRAUCHSANWEISUNG</b>	<b>16</b>
<b>[es]</b>	<b>INSTRUCCIONES DE USO</b>	<b>24</b>
<b>[fr]</b>	<b>MODE D'EMPLOI</b>	<b>34</b>
<b>[it]</b>	<b>ISTRUZIONI PER L'USO</b>	<b>43</b>
<b>[nl]</b>	<b>GEbruiksaanwijzingen</b>	<b>52</b>

**INDICATIONS FOR USE:**

A tool to transmit energy (provided by the RITA® 1500X RF Generator) for use in electrosurgery (ablation) in soft tissue, liver, and bone metastases.

**DESCRIPTION:**

The StarBurst® XLi-enhanced Electrosurgical Device consists of an insulated primary trocar with nine needles and four temperature sensors at the distal end. The four temperature sensors determine the average temperature and will be displayed in "TEMP 5" on the 1500X. A 25 cm length StarBurst® XLie SemiFlex device is also available that can be bent, if desired, in order to fit in a CT gantry. The device can bend to a radius of approximately 5 cm.

**WARNINGS:**


1. Do not attach anything (i.e., clamps, etc.) to the Device. This may damage the insulation, which could contribute to patient injury.
2. For the non-SemiFlex device, do not bend or kink the trocar or the needles. This may cause damage and result in a non-functional device.
3. For SemiFlex Devices, over bending of the trocar to a radius smaller than 5 cm beyond a 90° curvature and/or kinking the device can damage the trocar and cause patient injury.
4. Do not twist or exert high forces on the Device while it is deployed in the tissue. This may cause the needles to break and remain in the tissue.
5. Do not remove the Device without ensuring that the needles are fully retracted within the trocar.
6. Do not attach anything (i.e. clamp, etc.) to the tubing set as this may result in a kinked or blocked fluid path, which may result in an unpredictable lesion.
7. Patients with peripheral vascular deficiency are at increased risk of thermal injury from Dispersive Electrodes.
8. Patients with frail skin are at increased risk of skin damage from the adhesive on the Dispersive Electrodes.
9. Do not use metal introducers that do not have insulation. RF energy can be transmitted from the electrode through the un-insulated metal introducers to the patient, causing inadvertent burns.
10. If the device is being used in a laparoscopic procedure, care must be taken to avoid a gas embolism.


**PRECAUTIONS:**

1. Having RF power on at the same time as infusion using a method different from these instructions may alter the path of the electrical energy away from target tissues.
2. Keep the voltage/power as low as possible to achieve the desired end effect.
3. If retraction of the device becomes difficult, do not exert additional force. Infuse a small amount of saline solution and gently work the needles loose by alternating between a slightly retracted and a deployed position (while holding the main body with one hand and the trocar (at the point of insertion) with the other hand). The saline solution will generally soften the ablated tissue surrounding the electrodes. Once the ablated tissue softens, the arrays can be worked loose from the tissue and fully retracted.
4. If the device is being used in a laparoscopic procedure, activation of the device when not in contact with target tissue may cause capacitive coupling.
5. Reuse of single-use devices creates a potential risk of patient or user infections.
6. Contamination of the device may lead to injury, illness or death of the patient.
7. Reprocessing may compromise the integrity of the device and / or lead to device failure.

**EQUIPMENT LIST:**


- RITA® Model 1500X Generator with V8.12 or higher software.
- StarBurst® XLi-enhanced or StarBurst® XLie SemiFlex Electrosurgical Device
- RITA ThermoPad™ with adaptor cable (Dispersive Electrodes)
- IntelliFlow Pump
- Reusable Occlusion Bed
- Communication Cable
- Sterile 0.9% saline bag (250 mL required)

**XLi-enhanced**

**XLie SemiFlex****INSTRUCTIONS FOR USE:**

The following is the recommended procedure for operating the RITA® StarBurst® XLI-enhanced Device (Device) with micro-infusion.

1. Review and become familiar with the operation of the Model 1500X RF Generator (RF Generator) and the IntelliFlow Pump (Pump). Refer to the instructions for use accompanying each product.
2. Apply the Dispersive Electrodes according to the package instructions (see ThermoPad Instructions for Use).
3. Obtain a minimum of 250 mL of sterile 0.9% saline.
4. Inspect the Device packaging. If the sterile barrier is compromised, do not use.
5. Using sterile technique, open the packaging tray by peeling off the lid. Remove the Device.
6. Inspect the Device for any damage. Do not use a Device that has been damaged.
7. Inspect the Tubing Set prior to use, do not use the tubing set if there is any evidence of damage (e.g., kinked, cracked, leaking, etc.).
8. Connect the Main Cable to the RF Generator. Note: The Main Cable connectors are keyed to match the mating connectors. Therefore, attaching the Main Cable to the RF Generator requires minimal force. If force is required, you may be incorrectly connecting the cable and could damage the pins on the connectors. Verify that the number of pins on the connectors matches and verify that the pins of the connectors are not bent before continuing.
9. Verify the connections between Tubing Set and Device are secure.
10. Load the Tubing Set onto the Reusable Occlusion Bed. Turn Pump on. See IntelliFlow Pump Instructions for Use for more details.

**Tubing Set installed on the Reusable Occlusion Bed**

11. Load the Tubing Set and the Reusable Occlusion Bed onto the pump head.

---

**Note: Do not pull the pump lever forward at this time.**

---

**Tubing Set and Reusable Occlusion Bed installed in the IntelliFlow Pump**

12. Verify tubing orientation and flow direction, refer to label on top of Reusable Occlusion Bed.

**Note:** The pump lever will not lock into place if the tubing set has been loaded incorrectly.


13. Remove cap on Tubing Set spike and insert spike in the saline bag.

**Note:** Be sure to use sterile 0.9% normal injectable saline solution (at least 250 mL bag required).

14. Hang the saline bag on the IV pole mounted on the IntelliFlow Pump.
15. Purge the system of all air bubbles by allowing fluid to free flow through system to the pump.
16. Once the solution fills the Tubing Set, pull the pump lever forward securing the Reusable Occlusion Bed onto pump head.
17. For pump operation and use of the entire system refer to the IntelliFlow Pump instructions for use and the User's Guide and Service Manual for the 1500X RF Generator.
18. Verify the Generator is in the Purge Mode.
19. Press A on the Generator to purge and verify continuous flow of fluid from four of the microinfusion needles as well as the end of the trocar, press B to stop purge.
20. Verify that each of the Device temperature needles are functioning by holding the temperature needles between your sterile, gloved fingers. The temperature readings on the RF Generator should increase. If they do not, check connections and try again. Refer to the following diagram for location of temperature sensors.


**Note:** If flow is not observed out of four of the needles and the end of the trocar or if the thermocouples do not respond properly, the functionality of the device may be compromised.


#### Temperature Sensor Locations

21. Press the control mode on the generator to verify the Generator is in the XLi-enhanced mode.
22. If using an introducer system, ensure that the insulated introducer is from AngioDynamics, Inc. Refer to the introducer system instructions for use for introducer placement instructions.
23. Fully retract the Device needles by moving the slide on the handle backward (towards the cable connector) until it stops.
24. Using imaging guidance (e.g., ultrasound, CT), place the Device into the tissue. The tip of the trocar should be placed proximal to the center of the target area using the Placement Guide below.

Using the 1-cm markings on the trocar can assist in placement of the device.


**Note:** RF Generator will beep when the Device contacts the patient, indicating that low-level RF power is being delivered.

25. Deploy the Device needles by moving the slide forward (towards the trocar) on the handle to the desired final deployment size.

4-cm, 5-cm, 6-cm, 7-cm deployment markings

Unlabeled deployment markings represent 2-cm and 3-cm deployment

**Deployment Markings on Handle**

26. Verify position and deployment of the device with imaging (e.g., ultrasound, CT).
27. Retract to 2 cm deployment.
28. Verify that the Dispersive Electrodes are attached to the Generator, the Main Cable is connected to the Generator, and the Tubing Set is attached to the saline bag and pump and properly primed.
29. Program the RF Generator to the desired settings per the appropriate protocol listed below, and start RF power by depressing the foot pedal or the RF ON/OFF button. RF power will begin to increase.

#### Ablation of Soft Tissue or Liver Protocol:

**For a 7 cm ablation:** Power is preset to 250 W & Target Temperature is 105°.

Step	Deploy to:	Duration*
#1	2 cm	Maintain until Target Temp is reached (a beep will sound) then deploy to 3 cm.
#2	3 cm	Maintain until Target Temp is reached then deploy to 4 cm.
#3	4 cm	Maintain until Target Temp is reached then deploy to 5 cm.
#4	5 cm	Maintain until Target Temp is reached then deploy to 6 cm.
#5	6 cm	Set timer to 12 minutes. Maintain until 6 minutes remaining, then deploy to 7 cm.
#6	7 cm	Verify timer is at least 6 minutes (increase if necessary). When the time reaches 0.0 minutes the ablation cycle is complete. The system will automatically shut off the RF energy and go into "Cool Down Mode" for 30 seconds. At the end of the cool down mode temperatures above 60°C signal a successful ablation. (See section 30).

**For a 6 cm ablation:** Power is preset to 250 W & Target Temperature is 105°.

Step	Deploy to:	Duration*
#1	2 cm	Maintain until Target Temp is reached (a beep will sound) then deploy to 3 cm.
#2	3 cm	Maintain until Target Temp is reached then deploy to 4 cm.
#3	4 cm	Maintain until Target Temp is reached then deploy to 5 cm.
#4	5 cm	Set timer to 12 minutes. Maintain until 6 minutes remaining, then deploy to 6 cm.
#5	6 cm	Verify timer is at least 6 minutes (increase if necessary). When the time reaches 0.0 minutes the ablation cycle is complete. The system will automatically shut off the RF energy and go into "Cool Down Mode" for 30 seconds. At the end of the cool down mode temperatures above 60°C signal a successful ablation. (See section 30).

**For a 5 cm ablation:** Power is preset to 250 W & Target Temperature is 105°.

Step	Deploy to:	Duration*
#1	2 cm	Maintain until Target Temp is reached (a beep will sound) then deploy to 3 cm.
#2	3 cm	Maintain until Target Temp is reached then deploy to 4 cm.
#3	4 cm	Set timer to 12 minutes. Maintain until 6 minutes remaining, then deploy to 5 cm.
#4	5 cm	Verify timer is at least 6 minutes (increase if necessary). When the time reaches 0.0 minutes the ablation cycle is complete. The system will automatically shut off the RF energy and go into "Cool Down Mode" for 30 seconds. At the end of the cool down mode temperatures above 60°C signal a successful ablation. (See section 30).

**For a 4 cm ablation:** Power is preset to 250 W & Target Temperature is 105°.

Step	Deploy to:	Duration*
#1	2 cm	Maintain until Target Temp is reached (a beep will sound) then deploy to 3 cm.
#2	3 cm	Set timer to 12.0 minutes. Maintain until 6 minutes remaining, then deploy to 4 cm.
#3	4 cm	Verify timer is at least 6 minutes (increase if necessary). When the time reaches 0.0 minutes the ablation cycle is complete. The system will automatically shut off the RF energy and go into "Cool Down Mode" for 30 seconds. At the end of the cool down mode temperatures above 60°C signal a successful ablation. (See section 30).

#### Ablation of Painful Metastases with Extensive Bone Destruction Protocol:

**For a 5 cm ablation:** Set Target Temperature at 100° and Power at 150 W.

Step	Deploy to:	Duration*
#1	3 cm	Set timer to 20 minutes. When timer reaches 15 minutes, deploy to 4 cm.
#2	4 cm	Verify timer is at least 15 minutes (increase if necessary). When timer reaches 10 minutes, deploy to 5 cm.
#3	5 cm	Verify timer is at least 10 minutes (increase if necessary). When the time reaches 0.0 minutes the ablation cycle is complete. The system will automatically shut off the RF energy and go into "Cool Down Mode" for 30 seconds. At the end of the cool down mode temperatures above 60°C signal a successful ablation. (See section 30).

**For a 4 cm ablation:** Set Target Temperature at 100° and Power at 150 W.

Step	Deploy to:	Duration*
#1	3 cm	Set timer to 12 minutes. When timer reaches 8 minutes, deploy to 4 cm.
#2	4 cm	Verify timer is at least 8 minutes (increase if necessary). When the time reaches 0.0 minutes the ablation cycle is complete. The system will automatically shut off the RF energy and go into "Cool Down Mode" for 30 seconds. At the end of the cool down mode temperatures above 60°C signal a successful ablation. (See section 30).

**For a 3 cm ablation:** Set Target Temperature at 100° and Power at 150 W.

Step	Deploy to:	Duration*
#1	3 cm	Set timer to 8 minutes. When the time reaches 0.0 minutes the ablation cycle is complete. The system will automatically shut off the RF energy and go into "Cool Down Mode" for 30 seconds. At the end of the cool down mode temperatures above 60°C signal a successful ablation. (See section 30).

**For a 2 cm ablation:** Set Target Temperature at 100° and Power at 150 W.

Step	Deploy to:	Duration*
#1	3 cm	Set timer to 3 minutes. When the time reaches 0.0 minutes the ablation cycle is complete. The system will automatically shut off the RF energy and go into "Cool Down Mode" for 30 seconds. At the end of the cool down mode temperatures above 60°C signal a successful ablation. (See section 30).

**Ablation of Painful Metastases Encased in/Surrounded by Bone Protocol:**
**For a 5 cm ablation:** Set Target Temperature at 100° and Power at 150 W.

Step	Deploy to:	Duration*
#1	3 cm	Set timer to 12 minutes. When timer reaches 9 minutes, deploy to 4 cm.
#2	4 cm	Verify timer is at least 9 minutes (increase if necessary). When timer reaches 6 minutes, deploy to 5 cm.
#3	5 cm	Verify timer is at least 6 minutes (increase if necessary). When the time reaches 0.0 minutes the ablation cycle is complete. The system will automatically shut off the RF energy and go into "Cool Down Mode" for 30 seconds. At the end of the cool down mode temperatures above 60°C signal a successful ablation. (See section 30).

**For a 4 cm ablation:** Set Target Temperature at 100° and Power at 150 W.

Step	Deploy to:	Duration*
#1	3 cm	Set timer to 7 minutes. When timer reaches 5 minutes, deploy to 4 cm.
#2	4 cm	Verify timer is at least 5 minutes (increase if necessary). When the time reaches 0.0 minutes the ablation cycle is complete. The system will automatically shut off the RF energy and go into "Cool Down Mode" for 30 seconds. At the end of the cool down mode temperatures above 60°C signal a successful ablation. (See section 30).

**For a 3 cm ablation:** Set Target Temperature at 100° and Power at 150 W.

Step	Deploy to:	Duration*
#1	3 cm	Set timer to 3 minutes. When the time reaches 0.0 minutes the ablation cycle is complete. The system will automatically shut off the RF energy and go into "Cool Down Mode" for 30 seconds. At the end of the cool down mode temperatures above 60°C signal a successful ablation. (See section 30).

**For a 2 cm ablation:** Set Target Temperature at 100° and Power at 150 W.

Step	Deploy to:	Duration*
#1	3 cm	Set timer to 0.1 minute. When the time reaches 0.0 minutes the ablation cycle is complete. The system will automatically shut off the RF energy and go into "Cool Down Mode" for 30 seconds. At the end of the cool down mode temperatures above 60°C signal a successful ablation. (See section 30).

**\* Note: If the average temperature falls below the target temperature by more than 5 degrees for more than 5 seconds, the timer will discontinue counting down. Once the target temperature is reached again, the timer will resume counting down.**

30. Once the generator timer counts down to 0.0 minutes, the generator automatically enters into the Cool Down mode for 30 seconds. Cool down is a measure of how robust the ablation is. Watch the temperatures during the Cool Down mode. If the average temperature of the ablation after 30 seconds is below 55°C then it is more likely that the ablation is under sized than if the cool down temperature is over 55°C. If one array is significantly below 55°C then that area of the ablation is likely to be under developed and depending on the clinical need, additional ablation in that area is recommended. If temperatures are below 55°C consider retreating the area for an additional 5 minutes. Alternatively if 1 temperature is significantly lower than the others retract the tines, rotate the device 45 degrees, redeploy the tines and treat for an additional 5 minutes.
31. Additionally, confirm that the 30 second cool down temperatures of Temperatures 1 through 3 are 60°C or greater. If not then turn on RF power for 5 more minutes.
32. Upon completion of the desired ablation, turn the RF power off.
33. Fully retract the Device needles by pushing the slide on the handle backward (towards the cable connector) until it stops.
34. If desired, track ablation using 25 - 50 W can be performed during removal of the trocar.

**For each additional ablation:**

- Verify the continuous flow of fluid from four of the microinfusion needles as well as the end of the trocar between each ablation.
- Verify the geometry of the array of needles between each ablation.
- Gently remove any accumulated tissue from the device with wet gauze.
- 35. After completion of procedure, remove Tubing Set from reusable occlusion bed and dispose of Tubing Set and the Device.

**TROUBLESHOOTING GUIDE:****• If a needle becomes blocked:**

Turn off the RF energy and check that the fluid tubing has not become kinked, pinched or clamped.

**AND** gently remove any accumulated tissue from the device by wiping off the individual arrays with a coarse sponge or gauze. Do not use sharp objects to clean the needles as this can damage them and make the device not function properly.

**AND** check that tubing set is properly loaded on reusable occlusion bed. Reload the tubing onto the occlusion bed if necessary.

**AND** check for loose connections.

**AND** make sure tubing is properly purged (i.e. there are no large bubbles or air gaps inside the tubing).

**• If one temp is very different from the others:**

**IF** one temperature reading is low and the temperature needle is known to be located in or near a large vessel or outside the organ or target location, consider removing the unique temperature from the algorithm.

**IF** one temp is very low, but efficiency is okay (6-10), then consider leaving the green DEVICE TEMPERATURES (°C) LED on.

**IF** one temp is very low and efficiency is very low (0-5), consider removing the low temp from the algorithm by turning the green DEVICE TEMPERATURES (°C) LED off.

**IF** one temperature reading is very high for a short time (< 1 minute), consider leaving it as it is.

**IF** one temp is very high for more than a minute and efficiency is okay (6-10), consider taking the green DEVICE TEMPERATURES (°C) LED on out of the algorithm to allow the power to increase bringing the other temperatures up.

**• If efficiency is low (0-5):**

**IF** efficiency is low at the beginning of the case, check the dispersive electrodes for proper placement and ensure that the Device is fully deployed to desired ablation size.

**IF** efficiency is low during the case, consider lowering the target temperature, taking out the lowest temperature, or retracting the array and rotating the device.

**IF** power is high, consider decreasing the power.

**• If the “INFUSING – POWER ON in: XX” message appears on the generator:**

Impedance is too high and the system is automatically increasing the infusion rate to decrease the impedance while the RF power is turned off. RF power will resume once the impedance is within range. If this message appears more than three times, consider following the troubleshooting steps listed below for an impede-out event.

**• If there is an impede out (impedance too high):**

**IF** there is an impede out at the start of the ablation, check all connections and restart.

**IF** there is an impede out in the middle of an ablation, and the EFFICIENCY was gradually decreasing, consider retracting the array, rotating 45 degrees, redeploying, and continuing the ablation.

**IF** there is an impede out in the middle of an ablation, and the EFFICIENCY decreases sharply, check all connections and consider rotating and continuing the ablation.

**IF** there is an impede out at the end of an ablation, check cool down temperatures to determine if continued ablation is necessary.

**• If cool down temps are low:**

Consider ablating for an additional 5 minutes.

**Returns:** Defective devices may be returned to AngioDynamics Inc.  
AngioDynamics Customer Service should be contacted for all returns.

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© 2010, AngioDynamics, Inc.
RITA, StarBurst, and AngioDynamics are registered trademarks of AngioDynamics, Inc.
ThermoPad is a trademark of AngioDynamics, Inc.