


NE-STAR
Triple-Lumen Catheter



Performance Under Pressure

The Neostar* triple-lumen catheter's design gives clinicians multiple treatment options, consistent flow rates, and performance under pressure.



Description	UPN	Units/Box
Neostar Triple Lumen Central Venous Silicone Catheter with Cuff	H787CV332EK0	1

KIT COMPONENTS: (1) Triple Lumen Central Venous Silicone Catheter with Cuff; (1) 18 Ga Introducer Needle, extra thin wall, 7 cm length; (1) Guidewire, 0.038" OD x 76 cm length; (1) Vessel Dilator, 10F; (1) Peel Pro* PTFE Introducer, 13F; (1) Tunneler; (3) Injection Caps; (1) Syringe, 10cc; (1) #11 Scalpel; (6) Gauze Sponges; (1) Wound Dressing

IMPORTANT RISK INFORMATION

The following is a brief summary of important risk information for the AngioDynamics central venous catheter. For detailed information on the categories referenced, please consult the instructions for use packaged with each device. Observe all instructions prior to use. Failure to do so may result in patient complications.

INDICATIONS FOR USE: All central venous catheters are intended for long term vascular access. They are suitable for chemotherapy, infusion of IV fluids or drugs, parenteral nutrition, and transfusion or sampling of blood products.

CONTRAINDICATIONS: The central venous catheter is intended for long term vascular access only and should not be used for any purpose other than indicated in the instructions for use. Percutaneous catheter insertion in presence of thrombosis vessels, chronic obstructive lung disease, severe infection, or past irradiation of the subclavian area carries an increased risk.

WARNINGS AND PRECAUTIONS: Please see package insert for complete list of warnings and precautions.

POTENTIAL COMPLICATIONS: Consult package insert for a complete list of potential complications.

CAUTION: Federal (USA) law restricts these devices to sale by or on the order of a physician.



USA > 14 Plaza Drive, Latham, NY 12110 > tel: 800-772-6446 or 518-798-1215 > fax: 518-798-1360
 International > Haaksbergweg 75 (Margrietoren), 1101 BR, Amsterdam Z-O > The Netherlands
 tel: +31 (0)20 753 2949 > fax: +31 (0)20 753 2939

www.angiodynamics.com

*AngioDynamics, the AngioDynamics logo, NeoStar and Peel Pro are trademarks and/or registered trademarks of AngioDynamics, Inc., an affiliate or a subsidiary.

© 2014 AngioDynamics, Inc. MLC 499 US Rev B