

Designed by Clinicians for Clinicians



40ml of Saline
New Safety Cover for
Easy Vial Break

Safety Sharps Container



New Layout for
Increased Nursing
Efficiency

2 Chlorapreps
Tinted and Clear for
Site Preparation

Premium Safety Scalpel
and Needle



Extra Large Drapes
Ultra Absorbent Material

MORPHEUS®
SMART PICC
PREMIUM BEDSIDE INSERTION KIT



ANGIODYNAMICS®

THE ANGIODYNAMICS MORPHEUS SMART PICC PREMIUM BEDSIDE INSERTION KIT INCLUDES:

SMART PICC Catheter	1 Each	Drape, Absorbent	1 Each
Masks	2 Each	Lidocaine HCL 1%, 5cc Ampule	1 Each
SafetyGlide® Needle, 25 Ga. X 1.6cm	1 Each	Wipe, Alcohol	1 Each
Gloves	1 Pair	Tape Measures, Paper	2 Each
Surgical Tape	1 Each	Tourniquet	1 Each
Sodium Chloride (Saline) 0.9%, 10cc Ampules	4 Each	Non-Echo Needle, 21 Ga. X 7cm	1 Each
Drape, Fenestrated	1 Each	Polylined Towels	3 Each
Sponges, Gauze 10cm x 10cm (4in. x 4in.)	20 Each	Full Drape, 178cmx 254cm (70in. x 100in.)	1 Each
MaxPlus® Positive Displacement Connectors 1 per catheter lumen	1/lumen	CSR Wrap	1 Each
Sponges, Gauze 5cm x 5cm (2in.x 2in.)	6 Each	ChloraPrep® Antimicrobial Solution <i>Hi-Lite Orange® Applicator (1) Clear Applicator (1)</i>	2 Each
Statlock® Stabilization Device	1 Each	Syringe, Luer Lock 5ml	1 Each
Scissors, Disposable	1 Each	Bouffant Hair Nets	2 Each
Filter Straw® Aspiration Devices	2 Each	Sharps Container	1 Each
Syringes, Luer Lock 10ml	3 Each	Dressing, Transparent	1 Each
Guidewire, 0.46mm (0.018in.) x 80cm	1 Each	Premium Safety Scalpel	1 Each
Guidewire or Pre-loaded Stylet	1 Each		
Splittable 7 cm Sheath	1 Each		

Description	Catalog No.
4F Single Lumen CT PICC with Nitinol 0.018" Wire	12110902
5F Dual Lumen CT PICC with Nitinol 0.018" Wire	12110904
5F Dual PICC with 0.018" SS Wire	12110911
6F Dual PICC with 0.018" SS Wire	12110912
3F Single Lumen RT CT PICC with Nitinol 0.018" Wire	12110914
4F Single Lumen RT CT PICC with Nitinol 0.018" Wire	12110915
5F Single Lumen RT CT PICC with Nitinol 0.018" Wire	12110916
5F Dual Lumen RT CT PICC with Nitinol 0.018" Wire	12110917
6F Dual Lumen RT CT PICC with Nitinol 0.018" Wire	12110918
5F Dual Lumen RT CT PICC with 0.014"	12110919
3F Single Lumen CT PICC with Stylet	12110920
4F Single Lumen CT PICC with Stylet	12110921
5F Dual Lumen CT PICC with Stylet	12110923
3F Single Lumen RT CT PICC with Stylet	12110926

Description	Catalog No.
4F Single Lumen Reverse Taper CT PICC with Stylet	12110927
5F Single Lumen Reverse Taper CT PICC with Stylet	12110928
5F Dual Lumen Reverse Taper CT PICC with Stylet	12110929
4Fsingle RT PICC with 0.018" SS Wire	12110932
5F Dual RT PICC with 0.018" SS Wire	12110934
5F Dual RT PICC (55 cm) with Nitinol .018" Wire NEW!	12110936
5F Dual RT PICC (55 cm) with Stylet NEW!	12110937
5F Dual PICC (55 cm) with Nitinol .018" Wire NEW!	12110938
5F Dual PICC (55 cm) with Stylet NEW!	12110939
6F Triple Tapered PICC with Nitinol .018" Wire	12110941
6F Triple Tapered PICC with Stylet	12110942

RT = Reverse Taper

PHYSICIAN'S SIGNATURE

IMPORTANT RISK INFORMATION

INDICATIONS FOR USE:

The AngioDynamics Morpheus® SMART PICC CT is indicated for short or long term peripheral access to the central venous system for intravenous therapy, power injections of contrast media, and allows for central venous pressure monitoring. For blood sampling, infusion, or therapy use a 4 French or larger catheter. For central venous pressure monitoring, it is recommended that a catheter lumen of 20 gauge or larger be used.

CONTRAINDICATIONS:

This device is contraindicated whenever the patient

has or is suspected of having: infection, bacteremia, or septicemia; allergic reactions to the device's materials; insufficient body size to accommodate the size of the implanted device; past irradiation; previous episodes of venous thrombosis or vascular surgical procedures at the prospective placement site; local tissue factors will prevent proper device stabilization and/or access.

WARNINGS AND PRECAUTIONS:

Please see package insert for complete list of warnings and precautions.

POTENTIAL COMPLICATIONS:

Please see package insert for complete list of potential complications.

Indications, contraindications, warnings and instructions for use can be found in the instructions for use supplied with each device. Observe all instructions prior to use. Failure to do so may result in patient complications.

CAUTION: Federal (USA) law restricts these devices to sale by or on the order of a physician.

ANGIODYNAMICS®

USA > 603 Queensbury Avenue, Queensbury, NY 12804 > tel: 800-772-6446 or 518-798-1215 > fax: 518-798-1360
International > Building 2000, Beach Drive, IQ Cambridge, Waterbeach, Cambridge, CB25 9TE > United Kingdom
tel: +44 1223 729 300 > fax: +44 1223 729 349

www.angiodynamics.com

AngioDynamics and Morpheus are registered trademarks of AngioDynamics. Filter Straw is a registered trademark of B. Braun Medical Inc. Statlock is a registered trademark of C.R. Bard, Inc. or an affiliate. ChloraPrep and Hi-Lite Orange are registered trademarks of CareFusion 2200, Inc. SafetyGlide is a registered trademark of Becton, Dickinson & Co. MaxPlus is a registered trademark of Medegen, Inc. SPro is a trademark of InVaMed Technologies, Inc.

© 2011 AngioDynamics, Inc. MLC 233 Rev D 02/11