


DuraMax[®]

CHRONIC DIALYSIS CATHETER



Engineered for Excellence™

- Optimal Ease of Insertion
- High Flow–High Performance




Curved Tip
Catheter
TECHNOLOGY

ANGIODYNAMICS[®]

DuraMax® Engineered for Excellence™

CHRONIC DIALYSIS CATHETER



Description	Basic Kit		Catheter Only	
	Units/Box	Catalog No.	Units/Box	Catalog No.
20 cm Straight	□ 5	10302801	□ 5	10302901
22 cm Straight	□ 5	10302802	□ 5	10302902
24 cm Straight	□ 5	10302803	□ 5	10302903
28 cm Straight	□ 5	10302804	□ 5	10302904
32 cm Straight	□ 5	10302805	□ 5	10302905

Description	Basic Kit		Catheter Only	
	Units/Box	Catalog No.	Units/Box	Catalog No.
36 cm Straight	□ 5	10302806	□ 5	10302906
40 cm Straight	□ 5	10302807	□ 5	10302907
48 cm Straight	□ 5	10302808	□ 5	10302908
55 cm Straight	□ 5	10302809	□ 5	10302909

BASIC KIT INCLUDES: (1) DuraMax Dialysis Catheter, (1) 16F Valved Breakaway Introducer, (1) 12F Dilator, (1) 14F Dilator, (1) Tri-Ball Tunneler with Sleeve, (2) Injection Caps, (1) #11 Blade Safety Scalpel, (1) 0.038" J/Flex Guidewire, (1) Adhesive Wound Dressing, (1) 18 Ga. x 2-3/4" Introducer Needle

PHYSICIAN'S SIGNATURE

IMPORTANT RISK INFORMATION

INDICATIONS FOR USE: The DuraMax® Hemodialysis Catheter is indicated for use in attaining Long-Term vascular access for Hemodialysis and Apheresis. It may be inserted percutaneously and is primarily placed in the internal jugular vein of an adult patient. Alternate insertion sites include subclavian vein as required. Catheters greater than 40 cm are intended for femoral vein insertion. This catheter is indicated for > 30 days long term placement. The

Valved Peelable Introducer is intended for use in the percutaneous insertion of catheters in the venous system.

CONTRAINDICATIONS: The catheter is intended for Long-Term vascular access only and should not be used for any purpose other than indicated in these instructions. The valved peelable introducer sheath is not designed for use in the arterial system.

WARNINGS AND PRECAUTIONS: Please see package insert for complete list of warnings and precautions. **POTENTIAL COMPLICATIONS:** Consult package insert for a complete list of potential complications.

Indications, contraindications, warnings and instructions for use can be found in the instructions for use supplied with each device. Observe all instructions prior to use. Failure to do so may result in patient complications.

ANGIODYNAMICS®

USA > 603 Queensbury Avenue, Queensbury, NY 12804 > tel: 800-772-6446 or 518-798-1215 > fax: 518-798-1360
 International > Building 2000, Beach Drive, IQ Cambridge, Waterbeach, Cambridge, CB25 9TE > United Kingdom
 tel: +44 1223 729 300 > fax: +44 1223 729 349

www.angiodynamics.com

Covered under D603,044 with other US patents pending.
 AngioDynamics, DuraMax and Durathane are registered trademarks of AngioDynamics, Inc.
 Curved Tip and Engineered for Excellence are trademarks of AngioDynamics, Inc.
 SafeSheath D-Pro is a registered trademark of Pressure Products Medical Supplies, Inc.