


MORPHEUS®
SMART PICC

Ordering Information



ANGIODYNAMICS®

Insertion Kits

Catalog #	Catheter Description	Catheter Length	Guidewire
□ 12110902	4F Single Lumen Standard Taper	65 cm	.018 x 80 cm Nitinol Wire
□ 12110904	5F Dual Lumen Standard Taper	65 cm	.018 x 80 cm Nitinol Wire
□ 12110914	3F Single Lumen Reverse Taper	65 cm	.018 x 80 cm Nitinol Wire
□ 12110915	4F Single Lumen Reverse Taper	65 cm	.018 x 80 cm Nitinol Wire
□ 12110916	5F Single Lumen Reverse Taper	65 cm	.018 x 80 cm Nitinol Wire
□ 12110917	5F Dual Lumen Reverse Taper	65 cm	.018 x 80 cm Nitinol Wire
□ 12110918	6F Dual Lumen Reverse Taper	65 cm	.018 x 80 cm Nitinol Wire
□ 12110936	5F Dual Lumen Reverse Taper	55 cm	.018 x 80 cm Nitinol Wire
□ 12110938	5F Dual Lumen Standard Taper	55 cm	.018 x 80 cm Nitinol Wire
□ 12110941	6F Triple Lumen Reverse Taper	55 cm	.018 x 80 cm Nitinol Wire
□ 12110919	5F Dual Lumen Reverse Taper	65 cm	.014 x 80 cm Stainless Steel Wire
□ 12110911	5F Dual Lumen Standard Taper	65 cm	.018 x 80 cm Stainless Steel Wire
□ 12110912	6F Dual Lumen Standard Taper	65 cm	.018 x 80 cm Stainless Steel Wire
□ 12110932	4F Single Lumen Reverse Taper	65 cm	.018 x 80 cm Stainless Steel Wire
□ 12110934	5F Dual Lumen Reverse Taper	65 cm	.018 x 80 cm Stainless Steel Wire
□ 12110937	5F Dual Lumen Reverse Taper	55 cm	Stylet
□ 12110939	5F Dual Lumen Standard Taper	55 cm	Stylet
□ 12110942	6F Triple Lumen Reverse Taper	55 cm	Stylet
□ 12110920	3F Single Lumen Standard Taper	65 cm	Stylet
□ 12110921	4F Single Lumen Standard Taper	65 cm	Stylet
□ 12110923	5F Single Lumen Standard Taper	65 cm	Stylet
□ 12110926	3F Single Lumen Reverse Taper	65 cm	Stylet
□ 12110927	4F Single Lumen Reverse Taper	65 cm	Stylet
□ 12110928	5F Single Lumen Reverse Taper	65 cm	Stylet
□ 12110929	5F Dual Lumen Reverse Taper	65 cm	Stylet

Morpheus® SMART PICC Insertion Kit

KIT INCLUDES: (1) Smart PICC catheter, (2) Masks, (1) SafetyGlide® Needle, (1) Surgical Tape, (4) Sodium Chloride (Saline) 0.9%, 10cc Ampules, (1) Fenestrated Drape, (20) Sponges, Gauze 4 cm x 4cm, (1/lumen) MaxPlus® Positive Displacement Connectors, (6) Sponges, Gauze 10 cm x 10 cm, (1) Statlock® Stabilization Device, (1) Scissors, Disposable, (2) Filter Straw® Aspiration Devices, (3) Syringes, Luer Lock 10ml, (1) Guidewire, 0.46mm (0.018 in.) x 80 cm Guidewire or Pre-loaded Stylet, (1) Splittable 7cm Sheath, (1) Absorbent Drape, (1) Lidocaine HCL 1%, 5cc Ampule, (1) Alcohol Wipe, (2) Paper Tape Measures, (1) Tourniquet, (1) Non-Echo Needle, 21 Ga. X 7 cm, (3) Polylined Towels, (1) Full Drape, 178cm x 254cm (70in. x 100 in.), (1) CSR Wrap, (2) Chloraprep® Antimicrobial Solution – (1) Hi-Lite Orange® Applicator, (1) Clear Applicator, (1) Syringe, Luer Lock 5ml, (2) Bouffant Hair Nets, (1) Sharps Container, (1) Transparent Dressing, (1) Premium Safety Scalpel.



Standard Taper IR Kit

	Catheter Description	Stainless Steel Wire w/ Stainless Steel Tip	Nitinol Wire w/ Tungsten Tip	Stainless Steel Wire w/ Platinum Tip	
FULL IR KIT	3F Single Lumen CT PICC / .018 x 80 cm GW	□ 12100801	□ 12102801	□ 12104801	
	4F Single Lumen CT PICC / .018 x 80 cm GW	□ 12100802	□ 12102802	□ 12104802	
	5F Single Lumen CT PICC / .018 x 80 cm GW	□ 12100803	□ 12102803	□ 12104803	
	5F Dual Lumen CT PICC / .018 x 80 cm GW	□ 12100805	□ 12102805	□ 12104805	
	6F Dual Lumen CT PICC / .018 x 80 cm GW	□ 12100806	□ 12102806	□ 12104806	
	7F Single Lumen CT PICC / .018 x 130 cm GW	□ 12100807	□ 12102807	□ 12104807	
	3F Single Lumen CT PICC / .018 x 130 cm GW	□ 12100809	□ 12102809	□ 12104809	
	4F Single Lumen CT PICC / .018 x 130 cm GW	□ 12100810	□ 12102810	□ 12104810	
	5F Single Lumen CT PICC / .018 x 130 cm GW	□ 12100811	□ 12102811	□ 12104811	
	5F Dual Lumen CT PICC / .018 x 130 cm GW	□ 12100813	□ 12102813	□ 12104813	
	6F Dual Lumen CT PICC / .014 x 80 cm GW	□ 12100814	□ 12102814	□ 12104814	
	7F Dual Lumen CT PICC / .014 x 130 cm GW	□ 12100815	□ 12102815	□ 12104815	
	5F Dual Lumen CT PICC / .014 x 80 cm GW	□ 12100804	n/a	□ 12102804	
	5F Dual Lumen CT PICC / .014 x 130 cm GW	□ 12100812	n/a	□ 12102812	
	*5F Dual Lumen CT PICC / .018 x 80 cm GW	n/a	□ 12102835	n/a	
	*5F Dual Lumen CT PICC / .018 x 130 cm GW	n/a	□ 12102836	n/a	
	BASIC IR KIT	3F Single Lumen CT PICC / .018 x 80 cm GW	□ 12100817	□ 12102817	□ 12104817
		4F Single Lumen CT PICC / .018 x 80 cm GW	□ 12100818	□ 12102818	□ 12104818
5F Single Lumen CT PICC / .018 x 80 cm GW		□ 12100819	□ 12102819	□ 12104819	
5F Dual Lumen CT PICC / .018 x 80 cm GW		□ 12100821	□ 12102821	□ 12104821	
6F Dual Lumen CT PICC / .018 x 80 cm GW		□ 12100822	□ 12102822	□ 12104822	
7F Dual Lumen CT PICC / .018 x 80 cm GW		□ 12100823	□ 12102823	□ 12104823	
3F Single Lumen CT PICC / .018 x 130 cm GW		□ 12100825	□ 12102825	□ 12104825	
4F Single Lumen CT PICC / .018 x 130 cm GW		□ 12100826	□ 12102826	□ 12104826	
5F Single Lumen CT PICC / .018 x 130 cm GW		□ 12100827	□ 12102827	□ 12104827	
5F Dual Lumen CT PICC / .018 x 130 cm GW		□ 12100829	□ 12102829	□ 12104829	
6F Dual Lumen CT PICC / .018 x 130 cm GW		□ 12100830	□ 12102830	□ 12104830	
7F Dual Lumen CT PICC / .018 x 130 cm GW		□ 12100831	□ 12102831	□ 12104831	
5F Dual Lumen CT PICC / .014 x 80 cm GW		□ 12100820	n/a	□ 12102820	
5F Dual Lumen CT PICC / .014 x 130 cm GW		□ 12100828	n/a	□ 12102828	

*55 cm length catheter

Morpheus® SMART PICC Full IR Kit

Available in Standard and Reverse Taper PICC options.

FULL KIT INCLUDES: All Basic Kit components plus: one STATLOCK® stabilization device; one scalpel; one tape measure; one 10 mL syringe.

Morpheus® SMART PICC Basic IR Kit

Available with Standard PICC only.

BASIC KIT INCLUDES: One PICC; one guidewire; one sheath introducer; one 21ga. needle; two MaxPlus® positive displacement valves for dual lumen catheters; one MaxPlus® positive displacement valve for single lumen catheters.

Reverse Taper IR Kit

Catheter Description	Nitinol Wire w/ Tungsten Tip	Stainless Steel Wire w/ Platinum Tip
3F Single Lumen CT PICC / .018 x 80 cm GW	<input type="checkbox"/> 12102601	<input type="checkbox"/> 12104601
4F Single Lumen CT PICC / .018 x 80 cm GW	<input type="checkbox"/> 12102602	<input type="checkbox"/> 12104602
5F Single Lumen CT PICC / .018 x 80 cm GW	<input type="checkbox"/> 12102603	<input type="checkbox"/> 12104603
5F Dual Lumen CT PICC / .018 x 80 cm GW	<input type="checkbox"/> 12102605	<input type="checkbox"/> 12104605
6F Dual Lumen CT PICC / .018 x 80 cm GW	<input type="checkbox"/> 12102606	<input type="checkbox"/> 12104606
3F Single Lumen CT PICC / .018 x 130 cm GW	<input type="checkbox"/> 12102608	<input type="checkbox"/> 12104608
4F Single Lumen CT PICC / .018 x 130 cm GW	<input type="checkbox"/> 12102609	<input type="checkbox"/> 12104609
5F Single Lumen CT PICC / .018 x 130 cm GW	<input type="checkbox"/> 12102610	<input type="checkbox"/> 12104610
5F Dual Lumen CT PICC / .018 x 130 cm GW	<input type="checkbox"/> 12102612	<input type="checkbox"/> 12104612
6F Dual Lumen CT PICC / .018 x 130 cm GW	<input type="checkbox"/> 12102613	<input type="checkbox"/> 12104613
5F Dual Lumen CT PICC / .014 x 80 cm GW	n/a	<input type="checkbox"/> 12102604
5F Dual Lumen CT PICC / .014 x 130 cm GW	n/a	<input type="checkbox"/> 12102611
6F Triple Lumen CT PICC / .018 x 80 cm GW	<input type="checkbox"/> 12102617	n/a
6F Triple Lumen CT PICC / .018 x 130 cm GW	<input type="checkbox"/> 12102618	n/a
*5F Dual Lumen CT PICC / .018 x 80 cm GW	<input type="checkbox"/> 12102615	n/a
*5F Dual Lumen CT PICC / .018 x 130 cm GW	<input type="checkbox"/> 12102616	n/a

*55 cm length catheter

Reverse Taper Kits are only available in Full Kit options.

PHYSICIAN'S SIGNATURE

IMPORTANT RISK INFORMATION

INDICATIONS FOR USE:

The AngioDynamics Morpheus® SMART PICC CT is indicated for short or long term peripheral access to the central venous system for intravenous therapy, power injections of contrast media, and allows for central venous pressure monitoring. For central venous pressure monitoring, it is recommended that a catheter lumen of 20 gauge or larger be used.

CONTRAINDICATIONS:

This device is contraindicated whenever the patient has or is suspected of having: infection, bacteremia,

or septicemia; allergic reactions to the device's materials; insufficient body size to accommodate the size of the implanted device; past irradiation; previous episodes of venous thrombosis or vascular surgical procedures at the prospective placement site; local tissue factors will prevent proper device stabilization and/or access.

WARNINGS AND PRECAUTIONS:

Please see package insert for complete list of warnings and precautions.

POTENTIAL COMPLICATIONS:

Please see package insert for complete list of potential complications.

Indications, contraindications, warnings and instructions for use can be found in the instructions for use supplied with each device. Observe all instructions prior to use. Failure to do so may result in patient complications.

CAUTION: Federal (USA) law restricts these devices to sale by or on the order of a physician.

ANGIODYNAMICS®

USA > 603 Queensbury Avenue, Queensbury, NY 12804 > tel: 800-772-6446 or 518-798-1215 > fax: 518-798-1360
International > Building 2000, Beach Drive, IQ Cambridge, Waterbeach, Cambridge, CB25 9TE > United Kingdom
tel: +44 1223 729 300 > fax: +44 1223 729 349

www.angiodynamics.com

AngioDynamics and Morpheus are registered trademarks of AngioDynamics, Inc. Statlock is a registered trademark of C.R. Bard. Smart Taper is a trademark of AngioDynamics, Inc.

© 2011 AngioDynamics, Inc. MLC 235 Rev C 04/11