

DuraFlow

PRE-CURVE CHRONIC DIALYSIS CATHETER



DuraFlow

PRE-CURVE CHRONIC DIALYSIS CATHETER

Engineered for Durability



15.5 F DuraFlow[®] catheter facilitates high flow rates at modest arterial pressures. Lumen design provides improved kink resistance and flexibility in challenging tunnel locations.

Durathane[®] catheter material affords strength, softness, and resistance to a wide range of commonly used site care agents.



EmboSafe[®] peel-away sheath provides safety and security for vascular access and greatly reduces the risk of air embolism.

Oval guidewire holes allow for a smoother wire transition

Robust polyester cuff material allows optimal tissue ingrowth and catheter securement.

Rotating suture wings ensure patient comfort and flexibility to position catheter.

Description	Basic Kit		Full Kit		Catheter Only	
	Units/Box	UPN	Units/Box	UPN	Units/Box	UPN
20 cm Pre-Curve	□ 5	H787103012135	—	—	□ 10	H787103017105
22 cm Pre-Curve	—	—	—	—	□ 10	H787103017125
24 cm Pre-Curve	□ 5	H787103012065	—	—	—	—
28 cm Pre-Curve	□ 5	H787103012075	□ 5	H787103011075	—	—
32 cm Pre-Curve	□ 5	H787103012085	□ 5	H787103011085	—	—
36 cm Pre-Curve	□ 5	H787103012095	□ 5	—	□ 10	H787103017095

Description	Full Kit	
	Units/Box	UPN
24 cm Straight	□ 5	H787103011015
28 cm Straight	□ 5	H787103011025
32 cm Straight	□ 5	H787103011035
36 cm Straight	□ 5	H787103011045
40 cm Straight	□ 5	H787103011055

FULL KIT INCLUDES: (1) pair surgical gloves size 7-1/2, (4) 4" x 4" (10 cm) gauze sponges, (1) fenestrated drape with clear adhesive, (1) 3-pack providone iodine swab sticks, (2) 10 cc disposable luer lock syringes, (1) Raulerson bulb assembly, (1) 16F Valved peelable sheath, (1) DuraFlow™ dialysis catheter, (1) providone iodine ointment, (2) 5 cc disposable luer lock syringes, (1) 25 ga x 5/8" (1.6 cm) needle, (1) 22 ga x 1-1/2" (3.8 cm) needle, (1) straight disposable steel hemostat, (1) 18 ga x 1-1/2" (3.8 cm) needle, (1) 0.038" x 70 cm J/Flex guidewire, (1) #11 blade disposable scalpel with handle, (1) adhesive wound dressing, (2) injection caps, (1) 18 ga x 2-3/4" (7 cm) introducer needle, (1) tri-ball tunneler w/sleeve, (1) 12F vessel dilator, (1) 14F vessel dilator, (1) 2-0 silk suture w/curved needle, (1) 5 cc ampule lidocaine HCl 1%.

BASIC KIT INCLUDES: (1) DuraFlow Dialysis Catheter, (1) 16F Valved Breakaway Introducer, (1) 12F Dilator, (1) 14F Dilator, (1) Tri-Ball Tunneler with Sleeve, (2) Injection Caps, (1) #11 Blade Safety Scalpel, (1) 0.038" J/Flex Guidewire, (1) Adhesive Wound Dressing, (1) 18 Ga. x 2-3/4" Introducer Needle

IMPORTANT RISK INFORMATION

INDICATIONS FOR USE: The DuraFlow Hemodialysis Catheter is indicated for use in attaining long-term vascular access for hemodialysis and apheresis. It may be inserted percutaneously and is primarily placed in the internal jugular vein of an adult patient. Alternate insertion sites include subclavian vein as required. Catheters greater than 40 cm are intended for femoral vein insertion. This catheter is indicated for > 30 days long term placement. The

valved peelable introducer is intended for use in the percutaneous insertion of catheters in the venous system.

CONTRAINDICATIONS: The catheter is intended for long-term vascular access only and should not be used for any purpose other than indicated in these instructions. The valved peelable introducer sheath is not designed for use in the arterial system.

WARNINGS AND PRECAUTIONS: Please see package insert for complete list of warnings and precautions.

POTENTIAL COMPLICATIONS: Consult package insert for a complete list of potential complications.

Indications, contraindications, warnings and instructions for use can be found in the instructions for use supplied with each device. Observe all instructions prior to use. Failure to do so may result in patient complications.



USA > 14 Plaza Drive, Latham, NY 12110 > tel: 800-772-6446 or 518-798-1215 > fax: 518-798-1360
 International > Haaksbergweg 75 (Margrietoren), 1101 BR, Amsterdam Z-O > The Netherlands
 tel: +31 (0)20 753 2949 > fax: +31 (0)20 753 2939

www.angiodynamics.com

*AngioDynamics, the AngioDynamics logo, DuraFlow, EmboSafe and Durathane are trademarks and/or registered trademarks of AngioDynamics, Inc.

© 2015 AngioDynamics, Inc. MLC 457 US Rev C 09/15