



# SCHON

CHRONIC HEMODIALYSIS CATHETER



# Single Insertion Efficiency, Dual Catheter Effectiveness



> Over-the-wire insertion reduces risk of air embolism

**Polyurethane catheter**  
improves strength and flow rates

Priming marks and LOCK RIGHT\* adapters allows custom sizing

Single-site insertion

**HEMOLOCK\* anchoring hub**  
reduces bleeding and kinking

## POLYURETHANE SCHON CATHETER BASIC KIT

Description	Units/Box	UPN
14 cm Catheter	5	H787108018025
16 cm Catheter	5	H787108018035
18 cm Catheter	5	H787108018045
20 cm Catheter	5	H787108018055
22 cm Catheter	5	H787108018065
24 cm Catheter	5	H787108018075
26 cm Catheter	5	H787108018085

**BASIC KIT INCLUDES:** (1) Schon\* catheter, (2) LOCK RIGHT adapter sets, (2) injection caps, (2) temporary clamps

## SCHON INSERTION AND MEASUREMENT TRAY

Description	Units/Box	UPN
Cutting tip trocar (w/o catheter)	5	H787108023015
Blunt tip trocar (w/o catheter)	5	H787108023025

**KIT CONTENTS:** (1) 18g introducer needle, (1) J/Flex guidewire, (1) tear away sheath introducer, (1) #11 blade scalpel, (2) trocars, (1) 2-0 suture, (1) adhesive wound dressing, (2) injection caps, (2) extension sets, (2) clamps

## SCHON INSERTION AND MEASUREMENT TRAY REPLACEMENT PARTS

Description	Units/Box	UPN
Venous Lock Right adapters	10	10802401
Arterial Lock Right adapters	10	10802402
15 cm blunt tip trocars (tunnelers)	5	10800601
15 cm cutting tip trocars (tunnelers)	5	10800602
20 cm blunt tip trocars (tunnelers)	5	10800603

PHYSICIAN'S SIGNATURE

### IMPORTANT RISK INFORMATION

The following is a brief summary of important risk information for the AngioDynamics Schon Chronic Hemodialysis Catheter. For detailed information on the categories referenced, please consult the instructions for use packaged with each device. Observe all instructions prior to use. Failure to do so may result in patient complications.

**INDICATIONS FOR USE:** The AngioDynamics Schon Chronic Hemodialysis Catheter is designed for long term

hemodialysis and apheresis. It may be inserted percutaneously and is ideally placed in the internal jugular vein of an adult patient. Alternate insertion sites include subclavian vein as required.

**CONTRAINDICATIONS:** AngioDynamics Schon Chronic Hemodialysis Catheter is intended for long term vascular access only and should not be used for any purpose other than indicated in the instructions for use. Do not implant catheter in thrombosis vessels.

**WARNINGS AND PRECAUTIONS:** Please see package insert for complete list of warnings and precautions.

**POTENTIAL COMPLICATIONS:** Consult package insert for a complete list of potential complications.

**CAUTION:** Federal (USA) law restricts these devices to sale by or on the order of a physician.



USA > 14 Plaza Drive, Latham, NY 12110 > tel: 800-772-6446 or 518-798-1215 > fax: 518-798-1360  
International > Haaksbergweg 75 (Margrietoren), 1101 BR, Amsterdam Z-O > The Netherlands  
tel: +31 (0)20 753 2949 > fax: +31 (0)20 753 2939

[www.angiodynamics.com](http://www.angiodynamics.com)

\*AngioDynamics, the AngioDynamics logo and HEMOLOCK are trademarks and/or registered trademarks of AngioDynamics, Inc., an affiliate or a subsidiary. Schon is a trademark of Dr. Donald Schon. LOCK RIGHT is a trademark and/or registered trademark of Medical Components, Inc.

© 2015 AngioDynamics, Inc. MLC 035 Rev K 04/15