

Solero

Microwave Tissue Ablation System
SYSTEM SPECIFICATIONS



SIMPLE | SPEEDY | SCALABLE



SOLERO MICROWAVE TISSUE ABLATION SYSTEM

The Solero* Microwave Tissue Ablation (MTA) System and Accessories are indicated for the ablation of soft tissue during open procedures. The Solero MTA System is not intended for cardiac use.

SYSTEM SPECIFICATIONS

- 2.45 GHz Operating Frequency
- Solid State Generator
- Powers up to 140 W
- Intuitive Touch Screen User Interface
- Integrated Peristaltic Pump
- Applicator Coolant Temperature Monitoring
- Reflected Energy Monitoring
- No Grounding/ Dispersive/ Neutral Electrode Required
- Single Applicator System
- Applicator Available in 14, 19, and 29 cm Lengths

SKU	Description
H78712740000US0	Solero MW Generator
H787700106001US0	Solero MW Applicator 14 cm
H787700106002US0	Solero MW Applicator 19 cm
H787700106003US0	Solero MW Applicator 29 cm



Please refer to the Solero Generator Operator's Manual and the Solero Applicator Directions For Use for complete instructions, warnings, and precautions.



USA > 14 Plaza Drive, Latham, NY 12110 > tel: 800-772-6446 or 518-798-1215 > fax: 518-798-1360
International > Haaksbergweg 75 (Margriettoeren), 1101 BR, Amsterdam Z-O > The Netherlands >
tel: +31 (0)20 753 2949 > fax: +31 (0)20 753 2939

www.angiodynamics.com

*AngioDynamics, the AngioDynamics logo, Solero and the Solero logo are trademarks and/or registered trademarks of AngioDynamics, Inc., an affiliate or a subsidiary. © 2017 AngioDynamics, Inc. ANGB 358 US Rev 01 04/17