




## Percutaneous Introducer Sets

AngioDynamics percutaneous introducer sets, also known as Arnold Aspiration Sheaths, are part of a superior product offering designed to meet clinicians' vascular access needs.



18 cm Dilator assists in a smooth transition for easy insertion

10 cm Sheath includes side port extension tubing with a removable 3-way stopcock Hemo valve

Color coded hubs match standard French size for easy identification

Description	UPN	Hub Color
6F 18 cm Dilator, 10 cm Sheath, 18 Ga. Introducer Needle, .035/45 cm J-tip Guidewire	H787181002015	Green
6F 18 cm Dilator, 10 cm Sheath	H787181002035	Green
7F 18 cm Dilator, 10 cm Sheath	H787181002045	Orange

10 per box

#### IMPORTANT RISK INFORMATION

The following is a brief summary of important risk information for the AngioDynamics percutaneous introducer set. For detailed information on the categories referenced, please consult the instructions for use packaged with each device. Observe all instructions prior to use. Failure to do so may result in patient complications.

**INDICATIONS FOR USE:** AngioDynamics introducer products contain the essential components for the introduction of temporary pacing leads, balloon catheters, closed-end catheters or angiographic catheters into the vasculature during catheterization procedures.

**CONTRAINDICATIONS:** Introducers are not intended for critical care use.

**WARNINGS AND PRECAUTIONS:** Please see package insert for complete list of warnings and precautions.  
**POTENTIAL COMPLICATIONS:** Consult package insert for a complete list of potential complications.  
**CAUTION:** Federal (USA) law restricts these devices to sale by or on the order of a physician.



USA > 14 Plaza Drive, Latham, NY 12110 > tel: 800-772-6446 or 518-798-1215 > fax: 518-798-1360  
 International > Haaksbergweg 75 (Margrietoren), 1101 BR, Amsterdam Z-O > The Netherlands  
 tel: +31 (0)20 753 2949 > fax: +31 (0)20 753 2939

[www.angiodynamics.com](http://www.angiodynamics.com)