


# SCHON XL

ACUTE HEMODIALYSIS CATHETER



# Schon XL

## Acute Hemodialysis Catheter



Description	Catheter Only		Catheter Basic Set		Catheter Full Set	
	Units/Box	UPN	Units/Box	UPN	Units/Box	UPN
12 cm Acute Catheter	n/a	n/a	5	H787108007045	n/a	n/a
15 cm Acute Catheter	10	H787108017015	5	H787108007015	5	H787108027015
20 cm Acute Catheter	10	H787108017025	5	H787108007025	5	H787108027025
24 cm Acute Catheter	10	H787108017035	5	H787108007035	5	H787108027035

**BASIC SET INCLUDES:** (1) catheter, (1) 18g x 2-3/4" (7cm) introducer needle, (1) .035" x 70cm J/Flex guidewire, (1) #11 blade disposable scalpel with handle, (1) 2-0 suture w/curved needle, (1) HEMO-CATH clip, (1) adhesive wound dressing, (1) 12F vessel dilator, (1) 14F vessel dilator, (2) injection caps

**FULL SET INCLUDES:** (1) catheter, (1) pair surgical gloves, (4) 4"x4" (10cm) gauze sponges, (1) fenestrated drape with clear adhesive, (2) 10 mL disposable luer lock syringes, (1) 15F Vascu-Sheath Introducer, (2) 5 mL disposable luer lock syringes, (1) 25g x 5/8" (1.6cm) needle, (1) 22g x 1-1/2" (3.8cm) needle, (1) straight disposable steel hemostat, (1) 5 mL ampule Lidocaine HCL 1%, (1) 18g x 1-1/2" (3.8cm) needle, (1) .035" x 70cm J/Flex guidewire, (1) #11 blade disposable scalpel with handle, (1) adhesive wound dressing, (2) injection caps, (1) 18g x 2-3/4" (7cm) introducer needle, (1) HEMO-CATH clip, (1) 12F vessel dilator, (1) 14F vessel dilator, (1) 2-0 silk suture w/curved needle, (1) 3-pack Povidone Iodine swab sticks

PHYSICIAN'S SIGNATURE

### IMPORTANT RISK INFORMATION

The following is a brief summary of important risk information for the AngioDynamics Schon XL Acute Hemodialysis Catheter. For detailed information on the categories referenced, please consult the instructions for use packaged with each device. Observe all instructions prior to use. Failure to do so may result in patient complications.

**INDICATIONS FOR USE:** The AngioDynamics Schon XL Acute Hemodialysis Catheters is designed for acute

hemodialysis and apheresis. It may be inserted percutaneously and is ideally placed in the internal jugular vein. Although this catheter may be inserted into the subclavian or femoral vein, the internal jugular is the preferred site.

**CONTRAINDICATIONS:** AngioDynamics Schon XL Acute Hemodialysis Catheter is not intended for any use other than indicated in the instructions for use. Do not implant catheter in thrombosis vessels.

**WARNINGS AND PRECAUTIONS:** Please see package insert for complete list of warnings and precautions.

**POTENTIAL COMPLICATIONS:** Consult package insert for a complete list of potential complications.

**CAUTION:** Federal (USA) law restricts these devices to sale by or on the order of a physician.



USA > 14 Plaza Drive, Latham, NY 12110 > tel: 800-772-6446 or 518-798-1215 > fax: 518-798-1360  
 International > Haaksbergweg 75 (Margrietoren), 1101 BR, Amsterdam Z-O > The Netherlands  
 tel: +31 (0)20 753 2949 > fax: +31 (0)20 753 2939

[www.angiodynamics.com](http://www.angiodynamics.com)

\*AngioDynamics and the AngioDynamics logo are trademarks and/or registered trademarks of AngioDynamics, Inc., an affiliate or a subsidiary. Schon XL, HEMO-CATH and Vascu-Sheath are trademarks and/or registered trademarks of Medical Components, Inc.

© 2014 AngioDynamics, Inc. MLC 129 Rev H