

# ***EmboSafe***

VALVED, SPLITTABLE SHEATH DILATOR



## Reduce the Risk of Air Embolus and Blood Loss

AngioDynamics' EmboSafe<sup>®</sup> valved, splittable sheathed dilator is a valved introducer available for dialysis catheters. This proprietary valve reduces the risk of blood loss and air intake while allowing for smooth catheter insertion. The EmboSafe dilator requires no change in insertion technique, features a low profile sheath handle, and provides reliable peel performance. The EmboSafe dilator is available in sterile packages of five for 14F, 15F, and 16F sizes.



Description	Units/Box	UPN
EmboSafe 14F	5	H787065076015
EmboSafe 15F	5	H787065076025
EmboSafe 16F	5	H787065076035

PHYSICIAN'S SIGNATURE

### IMPORTANT RISK INFORMATION

**INDICATIONS FOR USE:** The EmboSafe valved, peelable introducer sheath is intended for use in the percutaneous insertion of pacing leads or catheters in the venous system.

**CONTRAINDICATIONS:** The EmboSafe valved peelable introducer sheath is not designed for use in the arterial system or as a hemostatic device.

**CAUTION:** Federal (USA) law restricts these devices to sale by or on the order of a physician.

**WARNINGS AND PRECAUTIONS:** The EmboSafe valved peelable introducer sheath is designed to reduce blood loss and the risk of air intake but it is not a hemostasis valve.

Please see package insert for complete list of warnings and precautions.



USA > 14 Plaza Drive, Latham, NY 12110 > tel: 800-772-6446 or 518-798-1215 > fax: 518-798-1360  
International > Haaksbergweg 75 (Margrietoren), 1101 BR, Amsterdam Z-O > The Netherlands  
tel: +31 (0)20 753 2949 > fax: +31 (0)20 753 2939

[www.angiodynamics.com](http://www.angiodynamics.com)

\*AngioDynamics, the AngioDynamics logo and EmboSafe are trademarks and/or registered trademarks of AngioDynamics, Inc., an affiliate or a subsidiary.

© 2014 AngioDynamics, Inc. MLC 165 Rev F