



Conventional Port Systems

TitanPort™



**PROVEN
PERFORMANCE**

TitanPort™

BENEFITS:

- > Blue boot locking mechanism
- > MRI compatible
- > Latex-free

TITANPORT™ TRAY COMPONENTS:

- > (1) TitanPort™ Port System
- > (1) Silicone Catheter
- > (2) Locking mechanisms (detached models)
- > (1) 22 Ga Non-coring needle
- > (1) 18 Ga Introducer needle
- > (1) PeelPro™ PTFE Introducer
- > (1) Guidewire, 0.035" x 50 cm
- > (2) 10 cc syringes
- > (1) Vein pick
- > (1) Blunt needle (detached models)
- > (1) Tunneler

TITANPORT™ KIT COMPONENTS:

- > (1) TitanPort™ Port System
- > (1) Silicone Catheter
- > (1) Locking mechanism (detached models)
- > (1) 22 Ga Non-coring needle
- > (1) Vein pick

The TitanPort™ line offers titanium models with silicone catheters.



Single Titanium Port	French	Part #
Detached silicone catheter with 8FR Introducer	7.2	P4355K
Attached silicone catheter with 8FR Introducer	7.2	P4305K
Detached silicone catheter with 10FR Introducer	9.6	P4455K
Attached silicone catheter with 10FR Introducer	9.6	P4405K