


The Technology Behind NeverTouch®


Featuring: Gold-Tipped Fiber

- Ultrasonically Visible
- Treats up to 4cm Veins



NEW NeverTouch-FRS

NeverTouch Fiber technology is now compatible with AngioDynamics' new Delta Laser systems as well as lasers formerly manufactured by Diomed.




 VenaCure EVLT™


ANGIODYNAMICS®




39.2 mm Great Saphenous Vein before a VenaCure® NeverTouch® procedure featuring the Gold-Tipped Fiber.



Vein after a VenaCure® NeverTouch® procedure featuring the Gold-Tipped Fiber with no flow under duplex.



Highly visible under ultrasound for ultimate security in placement and patient safety.



NeverTouch® fiber tip creates intervascular steam bubbles, while minimizing micro/macro perforations of vessel wall, that may occur if a fiber were to make contact with the vein wall.

References

1. Proebstle, *et al.* Nonocclusion and early reopening of the Great Saphenous Vein after Endovenous Laser Treatment is Fluence Dependent, *Dermatol Surg* 2004;30:174-178.
2. Proebstle, *et al.* Thermal Damage of the Inner Vein Wall During Endovenous Laser Treatment: key role of energy absorption by Intravascular Blood, *Dermatol Surg* 2002; 28:596-600.
3. Proebstle, *et al.* Endovenous treatment of the greater saphenous vein with a 940nm diode laser: Thrombotic occlusion after endoluminal thermal damage by laser-generated steam bubbles, *J Vasc Surg* 2002; 35:729-736.

Ordering Information

Catalog #	Description
<input type="checkbox"/> 11402001	NeverTouch Kit 45 cm
<input type="checkbox"/> 11402002	NeverTouch Kit 65 cm
<input type="checkbox"/> 11403001	NeverTouch-FRS Kit 45 cm
<input type="checkbox"/> 11403002	NeverTouch-FRS Kit 65 cm

Physician's Signature



AngioDynamics®, Inc.
 603 Queensbury Ave., Queensbury, NY 12804, U.S.A.
 Phone: 518-798-1215 Fax: 518-798-1360
 Toll-Free U.S.A.: 1-800-772-6446
www.angiodynamics.com www.venacureevlt.com

AngioDynamics, NeverTouch, and VenaCure are registered trademarks of AngioDynamics, Inc.
 VenaCure EVLT is a trademark of AngioDynamics.
 Covered by or for use under one or more of the following U.S. Patents: 6,769,433; 6,752,803; 6,258,084; 7,396,355; other U.S. Patents and foreign counterparts pending.

CAUTION: Federal Law (USA) restricts this device to sale by or on the order of a physician.

© AngioDynamics® Inc. 2009 MLC 230 Rev C 02/09