




Central Venous Catheters

NeoStar™ Triple Lumen Catheter



in a fluid environment, **BALANCE** is
ESSENTIAL.

Neostar™ Triple Lumen Catheter


BENEFITS:

- > Reinforced lumen walls
- > Consistent flow rates
- > Performance under pressure
- > Triple-lumen design—gives clinicians multiple treatment options.
- > The Neostar™ also features an improved, lower-profile fiber cuff to maximize tissue ingrowth
- > Silicone extension legs resist kinking from extensive clamping
- > One-piece ergonomic hub means easy, secure handling, mirror-smooth transition points. Soft outer web surfaces and flexible catheter enhance patient comfort

NEOSTAR™ KIT COMPONENTS:

- > (1) Triple Lumen Central Venous Silicone Catheter with Cuff
- > (1) 18 Ga Introducer Needle, extra thin wall, 7 cm length
- > (1) Guidewire, 0.038" OD x 76 cm length
- > (1) Vessel Dilator, 10FR
- > (1) Peel Pro™ PTFE Introducer, 13FR
- > (1) Tunneler
- > (3) Injection Caps
- > (1) Syringe, 10cc
- > (1) #11 Scalpel
- > (6) Gauze Sponges
- > (1) Wound Dressing

The Neostar™ triple-lumen catheter's revolutionary lumen design resists distortion under pressure, resulting in higher flows and lower pressures.



Description	Part #	Quantity
Neostar™ Triple Lumen Central Venous Silicone Catheter with Cuff	CV-332EK	1
Extension Leg Repair Kit	ELRK-400	1