

Micro-Access Tearaway


Introducers and Components

Evolving our products for all your needs.

AngioDynamics offers a wide selection of Introducer sizes, wire lengths and materials so you can choose a kit that fits your needs. These components can be used with the Morpheus[®] CT PICC products.

Our Introducer Kits include:

- 1 Introducer Needle
- 1 Guidewire
- 1 PTFE Micro-Access Tearaway Introducer



The Introducer

Choose from five Introducer sizes and 2 lengths. The tip is designed for ease of insertion and for a smooth transition. The PTFE sheath material provides a smooth peeling process.

The Introducer Needle

AngioDynamics Introducer needles allow for easy, less traumatic procedures.

Designed for Success

The .018" wire tips are atraumatic, designed to deflect off the vessel wall to ensure successful access. Available in Stainless and Platinum tip options.

Ordering Information

Introducer Kits

Catalog #	Introducer	Guidewire	Introducer Needle	Units/Box
<input type="checkbox"/> 06509901	3 Fr x 7 cm	.018" x 80 cm SS Mandrel/SSTip	21 Gauge x 7 cm Non-Echogenic	10
<input type="checkbox"/> 06509902	3 Fr x 10 cm	.018" x 80 cm SS Mandrel/SSTip	21 Gauge x 7 cm Non-Echogenic	10
<input type="checkbox"/> 06509903	3 Fr x 7 cm	.018" x 130 cm SS Mandrel/SSTip	21 Gauge x 7 cm Non-Echogenic	10
<input type="checkbox"/> 06509904	3 Fr x 10 cm	.018" x 130 cm SS Mandrel/SSTip	21 Gauge x 7 cm Non-Echogenic	10
<input type="checkbox"/> 06509905	3 Fr x 7 cm	.018" x 80 cm SS Mandrel/PtTip	21 Gauge x 7 cm Non-Echogenic	10
<input type="checkbox"/> 06509906	3 Fr x 10 cm	.018" x 80 cm SS Mandrel/PtTip	21 Gauge x 7 cm Non-Echogenic	10
<input type="checkbox"/> 06509907	3 Fr x 7 cm	.018" x 130 cm SS Mandrel/PtTip	21 Gauge x 7 cm Non-Echogenic	10
<input type="checkbox"/> 06509908	3 Fr x 10 cm	.018" x 130 cm SS Mandrel/PtTip	21 Gauge x 7 cm Non-Echogenic	10
<input type="checkbox"/> 06509909	4 Fr x 7 cm	.018" x 80 cm SS Mandrel/SSTip	21 Gauge x 7 cm Non-Echogenic	10
<input type="checkbox"/> 06509910	4 Fr x 10 cm	.018" x 80 cm SS Mandrel/SSTip	21 Gauge x 7 cm Non-Echogenic	10
<input type="checkbox"/> 06509911	4 Fr x 7 cm	.018" x 130 cm SS Mandrel/SSTip	21 Gauge x 7 cm Non-Echogenic	10
<input type="checkbox"/> 06509912	4 Fr x 10 cm	.018" x 130 cm SS Mandrel/SSTip	21 Gauge x 7 cm Non-Echogenic	10
<input type="checkbox"/> 06509913	4 Fr x 7 cm	.018" x 80 cm SS Mandrel/PtTip	21 Gauge x 7 cm Non-Echogenic	10
<input type="checkbox"/> 06509914	4 Fr x 10 cm	.018" x 80 cm SS Mandrel/PtTip	21 Gauge x 7 cm Non-Echogenic	10
<input type="checkbox"/> 06509915	4 Fr x 7 cm	.018" x 130 cm SS Mandrel/PtTip	21 Gauge x 7 cm Non-Echogenic	10
<input type="checkbox"/> 06509916	4 Fr x 10 cm	.018" x 130 cm SS Mandrel/PtTip	21 Gauge x 7 cm Non-Echogenic	10
<input type="checkbox"/> 06509917	5 Fr x 7 cm	.018" x 80 cm SS Mandrel/SSTip	21 Gauge x 7 cm Non-Echogenic	10
<input type="checkbox"/> 06509918	5 Fr x 10 cm	.018" x 80 cm SS Mandrel/SSTip	21 Gauge x 7 cm Non-Echogenic	10
<input type="checkbox"/> 06509919	5 Fr x 7 cm	.018" x 130 cm SS Mandrel/SSTip	21 Gauge x 7 cm Non-Echogenic	10
<input type="checkbox"/> 06509920	5 Fr x 10 cm	.018" x 130 cm SS Mandrel/SSTip	21 Gauge x 7 cm Non-Echogenic	10
<input type="checkbox"/> 06509921	5 Fr x 7 cm	.018" x 80 cm SS Mandrel/PtTip	21 Gauge x 7 cm Non-Echogenic	10
<input type="checkbox"/> 06509922	5 Fr x 10 cm	.018" x 80 cm SS Mandrel/PtTip	21 Gauge x 7 cm Non-Echogenic	10
<input type="checkbox"/> 06509923	5 Fr x 7 cm	.018" x 130 cm SS Mandrel/PtTip	21 Gauge x 7 cm Non-Echogenic	10
<input type="checkbox"/> 06509924	5 Fr x 10 cm	.018" x 130 cm SS Mandrel/PtTip	21 Gauge x 7 cm Non-Echogenic	10
<input type="checkbox"/> 06509933	6 Fr x 7 cm	.018" x 80 cm SS Mandrel/SSTip	21 Gauge x 7 cm Non-Echogenic	10
<input type="checkbox"/> 06509934	6 Fr x 10 cm	.018" x 80 cm SS Mandrel/SSTip	21 Gauge x 7 cm Non-Echogenic	10
<input type="checkbox"/> 06509935	6 Fr x 7 cm	.018" x 130 cm SS Mandrel/SSTip	21 Gauge x 7 cm Non-Echogenic	10
<input type="checkbox"/> 06509936	6 Fr x 10 cm	.018" x 130 cm SS Mandrel/SSTip	21 Gauge x 7 cm Non-Echogenic	10
<input type="checkbox"/> 06509937	6 Fr x 7 cm	.018" x 80 cm SS Mandrel/PtTip	21 Gauge x 7 cm Non-Echogenic	10
<input type="checkbox"/> 06509938	6 Fr x 10 cm	.018" x 80 cm SS Mandrel/PtTip	21 Gauge x 7 cm Non-Echogenic	10
<input type="checkbox"/> 06509939	6 Fr x 7 cm	.018" x 130 cm SS Mandrel/PtTip	21 Gauge x 7 cm Non-Echogenic	10
<input type="checkbox"/> 06509940	6 Fr x 10 cm	.018" x 130 cm SS Mandrel/PtTip	21 Gauge x 7 cm Non-Echogenic	10
<input type="checkbox"/> 06509941	7 Fr x 7 cm	.018" x 80 cm SS Mandrel/SSTip	21 Gauge x 7 cm Non-Echogenic	10
<input type="checkbox"/> 06509942	7 Fr x 10 cm	.018" x 80 cm SS Mandrel/SSTip	21 Gauge x 7 cm Non-Echogenic	10
<input type="checkbox"/> 06509943	7 Fr x 7 cm	.018" x 130 cm SS Mandrel/SSTip	21 Gauge x 7 cm Non-Echogenic	10
<input type="checkbox"/> 06509944	7 Fr x 10 cm	.018" x 130 cm SS Mandrel/SSTip	21 Gauge x 7 cm Non-Echogenic	10
<input type="checkbox"/> 06509945	7 Fr x 7 cm	.018" x 80 cm SS Mandrel/PtTip	21 Gauge x 7 cm Non-Echogenic	10
<input type="checkbox"/> 06509946	7 Fr x 10 cm	.018" x 80 cm SS Mandrel/PtTip	21 Gauge x 7 cm Non-Echogenic	10
<input type="checkbox"/> 06509947	7 Fr x 7 cm	.018" x 130 cm SS Mandrel/PtTip	21 Gauge x 7 cm Non-Echogenic	10
<input type="checkbox"/> 06509948	7 Fr x 10 cm	.018" x 130 cm SS Mandrel/PtTip	21 Gauge x 7 cm Non-Echogenic	10

Introducer Sets

Catalog #	Description	Units/Box
<input type="checkbox"/> 06511001	3 Fr x 10 cm Introducer	10
<input type="checkbox"/> 06511002	4 Fr x 10 cm Introducer	10
<input type="checkbox"/> 06511003	5 Fr x 10 cm Introducer	10
<input type="checkbox"/> 06511004	6 Fr x 10 cm Introducer	10
<input type="checkbox"/> 06511005	7 Fr x 10 cm Introducer	10
<input type="checkbox"/> 06511006	3 Fr x 7 cm Introducer	10
<input type="checkbox"/> 06511007	4 Fr x 7 cm Introducer	10
<input type="checkbox"/> 06511008	5 Fr x 7 cm Introducer	10
<input type="checkbox"/> 06511009	6 Fr x 7 cm Introducer	10
<input type="checkbox"/> 06511010	7 Fr x 7 cm Introducer	10

Guidewires

Catalog #	Description	Units/Box
<input type="checkbox"/> 06507210	.018" x 80 cm SS Mandrel/SSTip	10
<input type="checkbox"/> 06507211	.018" x 130 cm SS Mandrel/SSTip	10
<input type="checkbox"/> 06507212	.018" x 80 cm SS Mandrel/PtTip	10
<input type="checkbox"/> 06507213	.018" x 130 cm SS Mandrel/PtTip	10

Introducer Needles

Catalog #	Description	Units/Box
<input type="checkbox"/> 06508601	21 Gauge x 7 cm Echogenic Tip Needle	10
<input type="checkbox"/> 06508602	21 Gauge x 7 cm Non-Echogenic Tip Needle	10

Clinician's Signature

AngioDynamics, Inc. • 603 Queensbury Ave. • Queensbury, NY 12804 • USA
 Phone: 518-798-1215 • Fax: 518-798-1360 • Toll-Free USA: 1-800-772-6446
 International: AngioDynamics UK Ltd • Building 2000, Beach Drive, IQ Cambridge
 Waterbeach, Cambridge CB25 9TE • United Kingdom
 Tel: +44 1223 729 300 • Fax: +44 1223 729 349
www.angiodynamics.com

AngioDynamics and Morpheus are registered trademarks of AngioDynamics, Inc.
 CAUTION: Federal Law (U.S.) restricts the sale of this device by or on the order of a physician.
 ©AngioDynamics®, Inc. 2009 All rights reserved. MLC 232 REV B 01/09

ANGIODYNAMICS®

INCORPORATED