


TRANSJUGULAR ACCESS SET

THE...

...FINE NEEDLE SYSTEM




- ◆ Fine Needle System
- ◆ Fluoroscopic Guidance
- ◆ Improved Ross Cannula
- ◆ CO₂ or Iodine Compatible
- ◆ All PTFE Catheters

ANGIODYNAMICS®

INCORPORATED

TRANSJUGULAR ACCESS SET

Complete Set with Matched Components



*8 French catheter not shown above

Each Set Includes:

- 12 French Catheter (40cm long) with Hemostasis Valve and Side Arm
- 8 French* Curved PTFE Catheter (52cm long)
- 5 French PTFE Catheter (58cm long)
- 14 Gauge Modified Blunt Ross Needle (44cm long) with Orientation Handle
- 21 Gauge Puncture Needle (68cm long)
- 1cc Luer Lock Syringe
- 20cc Luer Lock Syringe
- Hemostasis Y-Adapter
- Dual Check Valve

ORDERING INFORMATION

Catalog Number	Description	Units/Box
15700101	Transjugular Access Set	1

ANGIODYNAMICS, INC.

603 Queensbury Ave.,
Queensbury, New York 12804 U.S.A.
Toll-Free: (800)-772-6446
Phone: (518)-798-1215
Fax: (518)-798-1360