

PULSE*SPRAY®/PRO™ Infusion Catheter TECHNICAL DESIGN

**Increased Infusion Ports
(SPRAY PATTERN COMPARISON)**

A second advantage of the slits is recognized in the ability to put more slits in a catheter than side holes without impacting the integrity of the catheter wall. This can be accomplished because no catheter wall material is removed. The PRO™ Infusion Catheter is able to incorporate 4 times the number of fluid ports found in common side hole catheters.

This difference can be visualized by comparing the following illustrations of active spray patterns.

