

E N G I N E E R E D  
F O R D U R A B I L I T Y

DURA-FLOW™ CHRONIC HEMODIALYSIS CATHETERS



ANGIODYNAMICS®

INCORPORATED

# DURA-FLOW™ SETTING THE STANDARD FOR HIGH PERFORMANCE,



## THE HYDRO-TIP® DESIGN

Maximizes flow, even when the catheter is against the vessel wall.




## OVER-THE-WIRE INSERTION

Reduces the risk of an air embolus. The oval holes allow for a smoother wire transition.



## THE ROTATING WING

Allows flexibility to position the catheter and rotate the wing for suturing.



Count on the 15.5F Dura-Flow™ Chronic Hemodialysis Catheter from AngioDynamics to deliver the performance-tough features you demand and the care your patient deserves. Engineered for durability, it combines all the attributes that are basic to the very best dialysis treatment and patient care.

## HIGH-TECH MATERIAL FOR EASIER CARE AND SITE MAINTENANCE

The Dura-Flow™ catheter is made of Durathane™—an advanced composite that outperforms traditional silicone and polyurethane while offering the advantages of both. It's exceptionally strong, biocompatible, and resistant to alcohol *and* iodine, for easier care and site maintenance.

## HIGH-FLOW DEPENDABLE

In addition to being one of the most durable catheters on the market, the Dura-Flow™ is non-position dependent. The Hydro-Tip® is specially designed to maximize flow rates — even when positioned against the vessel wall, resulting in a better dialysis treatment.

# STRENGTH AND VALUE.

Easy-to-read printed priming volumes

A rotating suture wing for flexible, secure anchoring and tip positioning

Polyester cuff anchors catheter



Flow-enhancing Hydro-Tip® design

## ORDERING INFORMATION

### DURA-FLOW™ FULL KIT

<i>Catalog No.</i>	<i>Description</i>	<i>Size</i>	<i>Units/Box</i>
10301114	Dura-Flow™ Straight Full	20 cm	5
10301112	Dura-Flow™ Straight Full	22 cm	5
10301101	Dura-Flow™ Straight Full	24 cm	5
10301102	Dura-Flow™ Straight Full	28 cm	5
10301103	Dura-Flow™ Straight Full	32 cm	5
10301104	Dura-Flow™ Straight Full	36 cm	5
10301105	Dura-Flow™ Straight Full	40 cm	5
10301111	Dura-Flow™ Straight Full	55 cm	5
10301113	Dura-Flow™ Pre-Curve Full	20 cm	5
10301110	Dura-Flow™ Pre-Curve Full	22 cm	5
10301106	Dura-Flow™ Pre-Curve Full	24 cm	5
10301107	Dura-Flow™ Pre-Curve Full	28 cm	5
10301108	Dura-Flow™ Pre-Curve Full	32 cm	5
10301109	Dura-Flow™ Pre-Curve Full	36 cm	5

Full Kit includes: (1) pair surgical gloves\* size 7-1/2, (4) 4" x 4" (10 cm) gauze sponges, (1) fenestrated drape with clear adhesive, (1) 3-pack providone iodine swab sticks, (2) 10 cc disposable luer lock syringes, (1) Raulerson bulb assembly, (1) 16F Valved peelable sheath, (1) Dura-Flow™ hemodialysis catheter, (1) providone iodine ointment, (2) 5 cc disposable luer lock syringes, (1) 25 ga x 5/8" (1.6 cm) needle, (1) 22 ga x 1-1/2" (3.8 cm) needle, (1) straight disposable steel hemostat, (1) 18 ga x 1-1/2" (3.8 cm) needle, (1) 0.038" x 70 cm J/Flex guidewire, (1) #11 blade disposable scalpel with handle, (1) adhesive wound dressing, (2) injection caps, (1) 18 ga x 2-3/4" (7 cm) introducer needle, (1) tri-ball tunneler w/sleeve, (1) 12F vessel dilator, (1) 14F vessel dilator, (1) 2-0 silk suture w/curved needle, (1) 5 cc ampule lidocaine HCl 1%.

\* CAUTION: THIS PRODUCT CONTAINS NATURAL RUBBER LATEX, WHICH MAY CAUSE ALLERGIC REACTIONS

### DURA-FLOW™ CATHETER ONLY

<i>Catalog No.</i>	<i>Description</i>	<i>Size</i>	<i>Units/Box</i>
10301711	Dura-Flow™ Straight Catheter Only	20 cm	10
10301713	Dura-Flow™ Straight Catheter Only	22 cm	10
10301701	Dura-Flow™ Straight Catheter Only	24 cm	10
10301702	Dura-Flow™ Straight Catheter Only	28 cm	10
10301703	Dura-Flow™ Straight Catheter Only	32 cm	10
10301704	Dura-Flow™ Straight Catheter Only	36 cm	10
10301705	Dura-Flow™ Straight Catheter Only	40 cm	10
10301714	Dura-Flow™ Straight Catheter Only	55 cm	10
10301710	Dura-Flow™ Pre-Curve Catheter Only	20 cm	10
10301712	Dura-Flow™ Pre-Curve Catheter Only	22 cm	10
10301706	Dura-Flow™ Pre-Curve Catheter Only	24 cm	10
10301707	Dura-Flow™ Pre-Curve Catheter Only	28 cm	10
10301708	Dura-Flow™ Pre-Curve Catheter Only	32 cm	10
10301709	Dura-Flow™ Pre-Curve Catheter Only	36 cm	10

### DURA-FLOW™ BASIC KIT

<i>Catalog No.</i>	<i>Description</i>	<i>Size</i>	<i>Units/Box</i>
10301214	Dura-Flow™ Straight Basic	20 cm	5
10301212	Dura-Flow™ Straight Basic	22 cm	5
10301201	Dura-Flow™ Straight Basic	24 cm	5
10301202	Dura-Flow™ Straight Basic	28 cm	5
10301203	Dura-Flow™ Straight Basic	32 cm	5
10301204	Dura-Flow™ Straight Basic	36 cm	5
10301205	Dura-Flow™ Straight Basic	40 cm	5
10301211	Dura-Flow™ Straight Basic	55 cm	5
10301213	Dura-Flow™ Pre-Curve Basic	20 cm	5
10301210	Dura-Flow™ Pre-Curve Basic	22 cm	5
10301206	Dura-Flow™ Pre-Curve Basic	24 cm	5
10301207	Dura-Flow™ Pre-Curve Basic	28 cm	5
10301208	Dura-Flow™ Pre-Curve Basic	32 cm	5
10301209	Dura-Flow™ Pre-Curve Basic	36 cm	5

Basic Kit includes: (1) Dura-Flow™ hemodialysis catheter, (1) #11 blade disposable scalpel with handle, (1) 0.038" x 70 cm J/Flex guidewire, (1) adhesive wound dressing, (1) 16F Valved peelable sheath, (1) tri-ball tunneler w/sleeve, (2) injection caps, (1) 18 ga x 2-3/4" (7 cm) introducer needle, (1) Raulerson bulb assembly, (1) 12F vessel dilator, (1) 14F vessel dilator.

### DIALYSIS ACCESSORIES

<i>Catalog No.</i>	<i>Description/Size</i>	<i>Units/Box</i>
05500901	Marked Guidewire 3 mm J Fixed Core (.038" x 70 cm uncoated, depth markings every 5 cm)	10
05500902	Marked Guidewire 3 mm J Fixed Core (.035" x 70 cm uncoated, depth markings every 5 cm)	10
10301403	DuraFlow™ Tunneler	5
10301801	Hemodialysis Catheter Repair Kit	5

**ANGIODYNAMICS®**  
INCORPORATED

ANGIODYNAMICS® INC., 603 Queensbury Ave., Queensbury, NY 12804, U.S.A.  
Phone: 518-798-1215 Fax: 518-798-1360 Toll-Free U.S.A.: 1-800-772-6446  
www.angiodynamics.com  
Dura-Flow™ and Durathane™ are trademarks of ANGIODYNAMICS®, INC.  
Hydro-Tip® and ANGIODYNAMICS® are registered trademarks of ANGIODYNAMICS, INC.  
MLC 151 REV D 02/07