

# EVENMORE®

Chronic Hemodialysis Catheter

## Strength Endurance Performance



ANGIODYNAMICS®

# ANGIODYNAMICS EVENMORE<sup>®</sup> Chronic Hemodialysis Catheter

Count on the 14.5F EVENMORE Catheter from AngioDynamics to deliver even greater performance and the care your patient deserves. Engineered for durability, it combines all the attributes that are basic to the very best dialysis treatment and patient care.

The 14.5F EVENMORE Catheter is the way to optimize blood flow and minimize venous back pressure. The low profile lumen to wall ratio in our dual lumen catheter is the key to enhanced flow. The round catheter shaft permits easier insertion and less patient discomfort.

The EVENMORE catheter is made of Durathane<sup>®</sup> — an advanced composite that outperforms traditional silicone and polyurethane while offering the advantages of both. It's exceptionally strong, biocompatible, and resistant to alcohol and iodine, for easier care and site maintenance.





Less than 5% recirculation, due to strategically placed venous return ports.



Luminal design allows for high flow.

# Providing even more durability and site care options.



# Ordering Information

	Catalog #	Description	Qty/Box
<b>EVENMORE® Basic Kit</b>	<input type="checkbox"/> 10303501	EVENMORE 24cm Straight Basic Kit	1
	<input type="checkbox"/> 10303502	EVENMORE 28cm Straight Basic Kit	1
	<input type="checkbox"/> 10303503	EVENMORE 32cm Straight Basic Kit	1
	<input type="checkbox"/> 10303504	EVENMORE 36cm Straight Basic Kit	1
	<input type="checkbox"/> 10303505	EVENMORE 40cm Straight Basic Kit	1
	<input type="checkbox"/> 10303506	EVENMORE 55cm Straight Basic Kit	1
	<input type="checkbox"/> 10303507	EVENMORE 24cm Pre-Curve Basic Kit	1
	<input type="checkbox"/> 10303508	EVENMORE 28cm Pre-Curve Basic Kit	1
	<input type="checkbox"/> 10303509	EVENMORE 32cm Pre-Curve Basic Kit	1
	<input type="checkbox"/> 10303510	EVENMORE 36cm Pre-Curve Basic Kit	1
	<input type="checkbox"/> 10303511	EVENMORE 20cm Straight Basic Kit	1
	<input type="checkbox"/> 10303512	EVENMORE 22cm Straight Basic Kit	1
	<input type="checkbox"/> 10303513	EVENMORE 20cm Pre-Curve Basic Kit	1
	<input type="checkbox"/> 10303514	EVENMORE 22cm Pre-Curve Basic Kit	1
<b>EVENMORE Catheter Only</b>	<input type="checkbox"/> 10303601	EVENMORE 24cm Straight Catheter Only	5
	<input type="checkbox"/> 10303602	EVENMORE 28cm Straight Catheter Only	5
	<input type="checkbox"/> 10303603	EVENMORE 32cm Straight Catheter Only	5
	<input type="checkbox"/> 10303604	EVENMORE 36cm Straight Catheter Only	5
	<input type="checkbox"/> 10303605	EVENMORE 40cm Straight Catheter Only	5
	<input type="checkbox"/> 10303606	EVENMORE 55cm Straight Catheter Only	5
	<input type="checkbox"/> 10303607	EVENMORE 24cm Pre-Curve Catheter Only	5
	<input type="checkbox"/> 10303608	EVENMORE 28cm Pre-Curve Catheter Only	5
	<input type="checkbox"/> 10303609	EVENMORE 32cm Pre-Curve Catheter Only	5
	<input type="checkbox"/> 10303610	EVENMORE 36cm Pre-Curve Catheter Only	5
	<input type="checkbox"/> 10303611	EVENMORE 20cm Straight Catheter Only	5
	<input type="checkbox"/> 10303612	EVENMORE 22cm Straight Catheter Only	5
	<input type="checkbox"/> 10303613	EVENMORE 20cm Pre-Curve Catheter Only	5
	<input type="checkbox"/> 10303614	EVENMORE 22cm Pre-Curve Catheter Only	5
<b>Dialysis Accessories</b>	<input type="checkbox"/> 06507602	EMBOSAFE™ Sheath	5
	<input type="checkbox"/> 05500901	Marked Guidewire 3mm J Fixed Core .038" x 70cm uncoated Depth markings every 5cm	10
	<input type="checkbox"/> 05500902	Marked Guidewire 3mm J Fixed Core .035" x 70cm uncoated Depth markings every 5cm	10
	<input type="checkbox"/> 10301401	Blunt Tip Tunneler	5
	<input type="checkbox"/> 10301402	Cutting Tip Tunneler	5
	<input type="checkbox"/> 10301801	Hemodialysis Catheter Repair Kit	5

Basic Kit Contents: (1) catheter, (1) valved peelable sheath introducer, (2) dilators, (1) tunneling tool, (1) scalpel, (1) J/Flex guidewire, (1) adhesive wound dressing, (1) introducer needle, (2) injection caps

Physician's Signature

**ANGIODYNAMICS®**

INCORPORATED

AngioDynamics, Inc. • 603 Queensbury Ave. • Queensbury, NY 12804, USA.  
Phone: 518-798-1215 • Fax: 518-798-1360 • Toll-Free USA: 1-800-772-6446  
[www.angiodynamics.com](http://www.angiodynamics.com)

AngioDynamics, Durathane and EvenMore are registered trademarks of AngioDynamics, Inc.  
EmboSafe is a trademark of AngioDynamics, Inc.  
Covered by one or more of the following US Patents: 6,712,789 and 4,808,155, foreign counterparts pending.

© AngioDynamics®, Inc. All rights reserved.

MLC 185 REV D 10/08