

# NiT-VU™

High-Performance Micro Guidewire

Have Visibility, Have Control



**ANGIODYNAMICS®**

INCORPORATED

## Super-Elastic Nitinol Shaft Construction

Provides excellent torque control, flexibility, kink resistance, control and maneuverability during placement

## Tungsten Wound Coil Tip

Provides excellent visibility during placement under fluoroscopy

## Highly Lubricious Coating



Reduces surface friction during wire advancement and placement

## Smooth Parabolic Tip Taper

Facilitates efficient, atraumatic placement through tortuous anatomy

	Catalog Number	Diameter (in)	Guidewire Length (cm)	Tip Stiffness	Tip Length (cm)	Tip Shape
<input type="checkbox"/>	05800118	0.014	80	Stiff	2	Angled
<input type="checkbox"/>	05800117	0.014	80	Stiff	2	Straight
<input type="checkbox"/>	05800120	0.014	180	Stiff	2	Straight
<input type="checkbox"/>	05800122	0.014	300	Stiff	2	Angled
<input type="checkbox"/>	05800119	0.014	80	Standard	5	Angled
<input type="checkbox"/>	05800104	0.014	80	Standard	5	Straight
<input type="checkbox"/>	05800121	0.014	180	Standard	5	Angled
<input type="checkbox"/>	05800105	0.014	180	Standard	5	Straight
<input type="checkbox"/>	05800123	0.014	300	Standard	5	Angled
<input type="checkbox"/>	05800106	0.014	300	Standard	5	Straight
<input type="checkbox"/>	05800101	0.014	80	Intermediate	7	Angled
<input type="checkbox"/>	05800102	0.014	180	Intermediate	7	Angled
<input type="checkbox"/>	05800103	0.014	300	Intermediate	7	Angled
<input type="checkbox"/>	05800125	0.018	80	Stiff	2	Angled
<input type="checkbox"/>	05800124	0.018	80	Stiff	2	Straight
<input type="checkbox"/>	05800128	0.018	180	Stiff	2	Angled
<input type="checkbox"/>	05800127	0.018	180	Stiff	2	Straight
<input type="checkbox"/>	05800130	0.018	300	Stiff	2	Straight
<input type="checkbox"/>	05800113	0.018	80	Standard	5	Angled
<input type="checkbox"/>	05800114	0.018	80	Standard	5	Straight
<input type="checkbox"/>	05800126	0.018	110	Standard	5	Straight
<input type="checkbox"/>	05800129	0.018	180	Standard	5	Angled
<input type="checkbox"/>	05800115	0.018	180	Standard	5	Straight
<input type="checkbox"/>	05800131	0.018	300	Standard	5	Angled
<input type="checkbox"/>	05800116	0.018	300	Standard	5	Straight
<input type="checkbox"/>	05800107	0.018	80	Intermediate	7	Angled
<input type="checkbox"/>	05800108	0.018	80	Intermediate	7	Straight
<input type="checkbox"/>	05800109	0.018	110	Intermediate	7	Straight
<input type="checkbox"/>	05800111	0.018	180	Intermediate	7	Angled
<input type="checkbox"/>	05800110	0.018	180	Intermediate	7	Straight
<input type="checkbox"/>	05800112	0.018	300	Intermediate	7	Angled

NOTES: Torque device included on all wires. Guidewires are packaged three per box.



\_\_\_\_\_  
Physician's Signature