

Maximal Protection for Improved Patient Care

Designed by nurses, the BioFlo PICC Maximal Barrier Nursing Kit includes a comprehensive set of components for safe and reliable bedside PICC placement.



BioFlo PICC
Maximal Barrier Nursing Kit

BIOFLO* PICC MAXIMAL BARRIER NURSING KIT INCLUDES:

BioFlo PICC	1 Each	Skin Protectant Swabstick	1 Each
Gown — XL	1 Each	Statlock* Stabilization Device	1 Each
Saline — Pre-filled 10 mL Syringes	4 Each	Dressing — Transparent	1 Each
Transducer Probe Cover Set (4 in x 58 in)	1 Each	Towel — Absorbent (14 in x 25 in)	2 Each
Drape Full Body (70 in x 100 in)	1 Each	Tourniquet	2 Each
Drape Fenestrated	1 Each	Surgical Tape	1 Each
Drape Underarm	1 Each	CSR Wrap	2 Each
Port Protector	1/Lumen	Mask with Earloops	2 Each
Tuohy Borst with Side Port	1 Each	Hair Net — Bouffant	2 Each
Stiffening Stylet	1 Each	Gloves — Surgical	1 Each
Gauze Pads (4 in x 4 in)	10 Each	Ampule Breaker Safety	1 Each
Introducer Safety Needle 21 Ga. x 2.75 in	1 Each	Lidocaine HCL 1%, 5 mL Ampule	1 Each
Introducer Needle 21 Ga. x 2.75 in — Super Sharp	1 Each	Label: Lidocaine	1 Each
Introducer Needle 21 Ga. x 1.57 in — Super Sharp (3F Single only)	1 Each	Syringe 5 mL	1 Each
Sharps Container	1 Each	Filter Straw	1 Each
Guidewire .018 x 45 cm Nitinol	1 Each	Needle 25 Ga. x .625 in	1 Each
Safety Scalpel — #11	1 Each	Tape Measure	2 Each
Sheath/Dilator Introducer	1 Each	Scissors	1 Each
Chloraprep* Clear 3 mL	1 Each	Patient Checklist	1 Each
Chloraprep Tinted Orange 3 mL	1 Each	Stop Sign	1 Each

Maximal Barrier Nursing Kit

Description	UPN	Catalog No.	Description	UPN	Catalog No.
BioFlo PICC 3F Single Lumen	H965750231	75-023	BioFlo PICC w/ PASV* Valve Technology 3F Single Lumen	H965750161	75-016
BioFlo PICC 4F Single Lumen	H965750121	75-012	BioFlo PICC w/ PASV Valve Technology 4F Single Lumen	H965750171	75-017
BioFlo PICC 5F Single Lumen	H965750131	75-013	BioFlo PICC w/ PASV Valve Technology 5F Single Lumen	H965750181	75-018
BioFlo PICC 5F Dual Lumen	H965750141	75-014	BioFlo PICC w/ PASV Valve Technology 5F Dual Lumen	H965750191	75-019
BioFlo PICC 6F Dual Lumen	H965750151	75-015	BioFlo PICC w/ PASV Valve Technology 6F Dual Lumen	H965750201	75-020
BioFlo PICC 6F Triple Lumen	H965750241	75-024	BioFlo PICC w/ PASV Valve Technology 6F Triple Lumen	H965750211	75-021
			BioFlo PICC w/ PASV Valve Technology 6F Hybrid Triple Lumen	H965750221	75-022

IMPORTANT RISK INFORMATION

BIOFLO PICC WITH ENDEXO TECHNOLOGY
 BIOFLO PICC WITH ENDEXO AND PASV VALVE TECHNOLOGY
 BIOFLO HYBRID PICC WITH ENDEXO AND PASV VALVE TECHNOLOGY

INTENDED USE/INDICATIONS FOR USE: The BioFlo PICCs are indicated for short or long-term peripheral access to the central venous system for intravenous therapy, including but not limited to, the administration of fluids, medications and nutrients; the sampling of blood; and for power injection of contrast media. Non-valved version is also indicated for central venous pressure monitoring.

CONTRAINDICATIONS: Venous thrombosis in any portion of the vein to be catheterized. Conditions that impede venous return from the extremity such as paralysis or lymphedema after mastectomy. Orthopedic or neurological conditions affecting the extremity. Anticipation or presence of dialysis grafts or other intraluminal devices. Hypercoagulopathy unless considerations are made to place the patient on anticoagulation therapy. Pre-existing skin surface or subsurface infection at or near the proposed catheter insertion site. Anatomical distortion of the veins from surgery, injury or trauma. Inadequate antecubital

veins. Anatomical irregularities (structural or vascular) which may compromise catheter insertion or catheter care procedures.

Indications, contraindications, warnings, potential complications and instructions for use can be found in the instructions for use supplied with each kit. Observe all instructions prior to use. Failure to do so may result in patient complications.

CAUTION: Federal (USA) law restricts these devices to sale by or on the order of a physician.



USA > 14 Plaza Drive, Latham, NY 12110 > tel: 800-772-6446 or 518-798-1215 > fax: 518-798-1360
 International > Haaksbergweg 75 (Margrietoren), 1101 BR, Amsterdam Z-O > The Netherlands
 tel: +31 (0)20 753 2949 > fax: +31 (0)20 753 2939

www.angiodynamics.com

Manufacturer:
 Navilyst Medical, Inc.,
 26 Forest Street, Marlborough, MA 01752

*AngioDynamics, the AngioDynamics logo, Navilyst Medical, BioFlo and PASV are trademarks and/or registered trademarks of AngioDynamics Inc., an affiliate or a subsidiary. Endexo is a trademark and/or registered trademark of Interface Biologics. Chloraprep is a trademark and/or registered trademark of Carefusion 2200, Inc. Statlock is a trademark and/or registered trademark of C. R. Bard, Inc.